INSTRUCTIONS FOR USE – MANUAL ENDODONTIC INSTRUMENTS

EN	Instructions for use	2
BG	Инструкции за употреба.....	7
CS	Návod k použití.....	12
DA	Brugsanvisning	17
DE	Gebrauchsanweisung	22
EL	Οδηγίες χρήσης	27
ES	Instrucciones de uso.....	32
ET	Kasutusjuhend.....	37
FI	Käyttöohjeet.....	42
FR	Notice d'utilisation	47
HR	Upute za uporabu.....	52
HU	Használati utasítás.....	57
IT	Istruzioni per l'uso	62
LT	Naudojimo instrukcija	67
LV	Lietošanas instrukcijas.....	72
NL	Gebruiksaanwijzing	77
NO	Bruksanvisning	82
PL	Instrukcja użytkowania	87
PT	Instruções de utilização.....	92
RO	Instrucţiuni de utilizare	97
SK	Návod na použitie	102
SL	Navodila za uporabo	107
SV	Bruksanvisning	112
TR	Kullanım talimatları	117

FOR DENTAL USE ONLY

INSTRUCTIONS FOR USE – MANUAL ENDODONTIC INSTRUMENTS

1) REFERENCES

References	Designation	Sizes	Sterile	Reusable*	Length in mm		
					21	25	31
S1.1B4.xx.xxx.FK	K-Files	006-010	YES	NO	X	X	X
		015-080	YES	YES	X	X	X
01.154.xx.xxx.CG	K-Files	006-010	NO	NO	X	X	X
		015-080	NO	YES	X	X	X
01.153.xx.xxx.CG	H-Files	008-010	NO	NO	X	X	X
		015-080	NO	YES	X	X	X
01.553.xx.xxx.CG	H-Files Flex	008-010	NO	NO	X	X	X
		015-080	NO	YES	X	X	X
01.160.xx.xxx.CG	Reamers	006-010	NO	NO	X	X	X
		015-080	NO	YES	X	X	X
03.050.31.xxx.CT	Barbed Broaches	025-060	NO	NO	X		

* Reusable devices bear the following 8 petaled symbol () at the distal end of the handle. These devices can be used and reprocessed up to 8 times, each petal corresponding to 1 canal.

2) PARAMETERS FOR USE

Manual endodontic instruments are intended for hand-held use only.

3) INTENDED USE

Manual endodontic instruments intended for root canal treatment.

4) INDICATIONS FOR USE

Manual endodontic instruments are indicated in case of pulp or periapical diseases that require endodontic treatment. With the exception of Barbed Broaches, instruments can be also used in case of retreatment (i.e. secondary treatment).

All instruments are designed for preparation and shaping of root canals.

Manual endodontic instruments are intended for use in medical or hospital facilities, by qualified health professionals.


5) KEY PERFORMANCES

Manual instruments are compliant with all applicable testing requirements of the ISO 3630 series of standards.

6) COMPOSITION

The working part of these instruments is made of Stainless Steel and carries a plastic handle. The manual endodontic instruments carry a silicone endo stop to determine the working length of root canal preparation, except for Barbed Broaches products. In addition, except for Barbed Broaches products, all instruments carry a safety hole.


Product example: K-File

Identification of instruments with coding according to the handle color code, symbol and tip size.
Identification of the working part length with silicone endo stop colors:

Product example: Barbed Broach

7) CONTRAINDICATION

No contraindication based on the patient population is known.

8) ADVERSE REACTIONS

None known.

9) WARNINGS AND PRECAUTIONS

For your safety, the safety of your patients and the safety of third parties, observe the safety notes below to minimize the following risks to the manual endodontic instruments:

User qualification and patient population

- The instructions for use are a component of the product and must be read carefully prior to use and be accessible at all times.
- The manual endodontic instruments may only be used in accordance with the intended use, any other type of use is not permitted.
- These instruments have not been tested on children, pregnant nor breastfeeding women.

Infection / Toxic or allergic reaction hazards

- Manual endodontic instruments must not be used on individuals with a known allergic sensitivity to stainless steel and silicone.
- For your own safety, wear personal protective equipment (masks, gloves, and goggles) during the treatment and reprocessing of the instrument.
- Do not use if sterile packaging seal is broken, damaged or wet.
- Do not use the instruments after expiration date
- Do not use if non-sterile packaging is broken.
- Non-sterile instruments should be reprocessed before the first use and after each use.
- Single use instruments must not be retreated for subsequent use. This could comprise their integrity.
- For reusable instruments, use approved cleaning and disinfecting agents:
 - Washing: Neodisher® Mediclean Forte (0.5 % concentration)
 - Thermal disinfection: Neodisher® Mediklar Special (0.03 % concentration)
 in accordance with the reprocessing instructions on FKG's website.
- Reprocess the instrument appropriately before disposal.

Breakage / Deterioration / Inhalation or swallowing hazards

A damaged instrument or components could injure patients, users and third parties.

- Carefully read the label, marking on the packaging and this instruction for use to ensure the correct identification and use of device.
- Use a dental dam and a protective wire inside the safety hole when using manual endodontic instruments to avoid aspiration or ingestion by the patient.
- Take several radiographs from different angles to determine the anatomy of the root canals (length, width and curvature).
- Always inspect the instrument(s) before use, and discard (it) them if there are any visible defect(s).
- Cracks, deformations, signs of corrosion, loss in color or marking are signs that the device is no longer able to achieve the required performance level and should be discarded.
- Regularly clean flutes and inspect instrument during use, and discard if it shows any signs of distortion or wear.
- If manual endodontic instrument is not progressing easily, withdraw it from the canal, clean it, check its cutting edges, then irrigate the canal and repeat the process.
- Do not exceed a sterilization temperature of 135°C.

Device packaging

- If the packaging has been opened, damaged, become wet or the expiry date has passed, the sterile state of the instruments inside the packaging is not guaranteed. Discard the device(s) that are single use and consult the reprocessing instructions available on FKG's website for multiple use device(s).
- The packaging of the instruments should be carefully examined before opening (packaging integrity, no humidity, no visual degradation of instruments) to ensure that the packaging's integrity has not been compromised during transport and storage.

Warnings and precautions for the (re)processing of instruments

- Manual endodontic instruments that are provided sterile, may be reprocessed for multi-usage following reprocessing instructions available on FKG's website, except for single use instruments (recognizable with following symbol ② on the product packaging) which must be used for a single canal only.
- Manual endodontic instruments that are provided non-sterile, must be reprocessed before the first use and between each use, following reprocessing instructions available on FKG's website, except for single use instruments (recognizable with following symbol ② on the product packaging) which must be used for a single canal only.

Furthermore, it is the responsibility of the practitioner to always check its instruments before each use to identify possible signs of defects.

In case of accident

- All serious events occurring in relation to the product must be reported to the manufacturer and the competent authority according to local regulations.

10) PROTOCOL FOR USE

- 1) Create straight-line coronal and radicular access.
- 2) If necessary, use Barbed Broach to remove pulp tissue. *Barbed Broach must always enter the root canal gently.*
- 3) Manual endodontic instruments (i.e. K-Files, H-Files, H-Files Flex, Reamers) may have several possible uses during a root canal treatment. Indeed, the instruments can be used for probing, exploration, patency confirmation, working length determination (using a radiograph or, if required, an apex locator), recapitulation, debris extraction, glide path, retreatment and/or root canal shaping. Manual endodontic instruments can also be used in association to mechanized instruments. *RACE® EVO, R-Motion® or XP-endo® Rise are recommended solutions for mechanized root canal treatments.*
 - a) K-files and Reamers are suitable for root canal preparation in filing and reaming motion, respectively. These instruments are intended for enlarging canals. Rotations must not exceed 360°.
 - b) Hedström manual endodontic instruments (i.e. H-Files, H-Files Flex) are intended only for translational strokes inside the canal. Rotational working movements are not recommended because of the possibility of fracture.
- 4) Repeat the appropriate process until the canal is sufficiently enlarged to proceed to the next size instrument. *Do not proceed too quickly - it is preferable to use an instrument at its maximum dimensions before proceeding to the next instrument in the series.* Irrigate (i.e. NaOCl) frequently and abundantly throughout the process. The instrument must be cleaned repeatedly during use; never use an instrument packed with filings or debris.

Once all the canals have been shaped, proceed with the cleaning of the canals. *For an optimal result, the use of Finisher is recommended.*

Once the cleaning of the canals is completed, rinsed and dried; proceed to the next phase of treatment:

5) Obturate with gutta-percha and sealer. TotalFill® BC Points™ and TotalFill® BC Sealer™ are recommended.

11) CONDITIONS AND TERM OF STORAGE, SHELF-LIFE

The shelf-life is 5 years for sterile instruments. The use-by date is indicated on the packaging.

Sterility is preserved and guaranteed for the whole shelf life unless the packaging is broken, or the storage conditions are not respected.

No specific shelf-life is indicated for the non-sterile instruments.

12) REPROCESSING PROTOCOL

Reprocessing must be applied to manual instruments before the first use (only if devices are supplied non-sterile, see section 1) and before each reuse (only if devices are reusable, see section 1). Instruments must be reprocessed following reprocessing instructions available on FKG's website.

Due to the design of the instruments and/or materials used and in the absence of contrary indications in the labeling or instructions for use of the device, the number of uses is maximum 8 canals.

13) DISPOSAL

When a device reaches the end of its life, make sure that it is discarded in accordance with the applicable laws and regulations.

14) SYMBOLS

For explanation of symbols for IFUs and labels, please consult the document "Glossary of the symbols used by FKG Dentaire" available on FKG's website.

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