## INSTRUCTIONS FOR USE – PASTE FILLERS

|           |                                 |            |
|-----------|---------------------------------|------------|
| <b>EN</b> | Instructions for use .....      | <b>2</b>   |
| <b>BG</b> | Инструкции за употреба .....    | <b>6</b>   |
| <b>CS</b> | Návod k použití .....           | <b>11</b>  |
| <b>DA</b> | Brugsanvisning .....            | <b>16</b>  |
| <b>DE</b> | Gebrauchsanweisung .....        | <b>20</b>  |
| <b>EL</b> | Οδηγίες χρήσης .....            | <b>25</b>  |
| <b>ES</b> | Instrucciones de uso .....      | <b>30</b>  |
| <b>ET</b> | Kasutusjuhend .....             | <b>35</b>  |
| <b>FI</b> | Käyttöohjeet .....              | <b>39</b>  |
| <b>FR</b> | Notice d'utilisation .....      | <b>43</b>  |
| <b>HR</b> | Upute za uporabu .....          | <b>48</b>  |
| <b>HU</b> | Használati utasítás .....       | <b>52</b>  |
| <b>IT</b> | Istruzioni per l'uso .....      | <b>57</b>  |
| <b>LT</b> | Naudojimo instrukcija .....     | <b>61</b>  |
| <b>LV</b> | Lietošanas instrukcijas .....   | <b>66</b>  |
| <b>NL</b> | Gebruiksaanwijzing .....        | <b>70</b>  |
| <b>NO</b> | Bruksanvisning .....            | <b>75</b>  |
| <b>PL</b> | Instrukcja użytkowania .....    | <b>79</b>  |
| <b>PT</b> | Instruções de utilização .....  | <b>84</b>  |
| <b>RO</b> | Instrucţiuni de utilizare ..... | <b>89</b>  |
| <b>SK</b> | Návod na použitie .....         | <b>94</b>  |
| <b>SL</b> | Navodila za uporabo .....       | <b>99</b>  |
| <b>SV</b> | Bruksanvisning .....            | <b>103</b> |
| <b>TR</b> | Kullanım talimatları .....      | <b>107</b> |

FOR DENTAL USE ONLY

## INSTRUCTIONS FOR USE – PASTE FILLERS

### 1) REFERENCES

| References       | Designation          | Sizes   | Sterile | Reusable* | Length in mm |    |    |
|------------------|----------------------|---------|---------|-----------|--------------|----|----|
|                  |                      |         |         |           | 21           | 25 | 31 |
| 02.171.XX.XXX.GM | Lentulo              | 025-040 | NO      | YES       | X            | X  | X  |
| 02.190.XX.XXX.GM | Spring Paste Fillers | 025-040 | NO      | YES       | X            | X  | X  |
| 02.180.XX.XXX.GS | Sensipast            | 025-040 | NO      | YES       | X            | X  | X  |

\* All instruments are reusable and can be used up to 10 times.

### 2) PARAMETERS FOR USE AND COMPATIBLE DEVICES

Paste fillers instruments are intended for use with an endodontic motor compliant with ISO 1797.

#### Lentulo

Speed: 800 rpm  
Torque limit: 1.5 Ncm

#### Spring Paste Fillers

Speed: 800 rpm  
Torque limit: 1.5 Ncm

#### Sensipast

Speed: 800 rpm  
Torque limit: N/A

No torque limit is needed for Sensipast because Sensipast is a paste filler with an automatic safety clutch. Sensipast reacts before the restraints exceed the metal elasticity limits.

### 3) INTENDED USE

Paste fillers instruments are intended for conveying filling material or medicaments into a root canal system.

#### 4) INDICATIONS FOR USE

Paste fillers instruments are indicated in case of dental pulp and/or dental periapical diseases that require endodontic treatment. All instruments are designed to ensure efficient and effective endodontic obturation technique.

Paste fillers instruments are intended for use in medical or hospital facilities, by qualified health professionals.


#### 5) KEY PERFORMANCES

Paste fillers instruments are compliant with all applicable testing requirements of the ISO 3630 series of standards.


#### 6) COMPOSITION

The working part of these instruments is made of Stainless Steel and carries a shank of Maillechort (for Lentulo and Spring Paste Fillers only) or Stainless Steel (for Sensipast only). Tip size is reported on shank by a colored ring (for Lentulo and Spring Paste Fillers only) or a colored obturator (for Sensipast only).


*Product example: Lentulo*


Identification of the ISO diameter with silicone ring colors:

*Product example: Sensipast*

Identification of the ISO diameter with silicone obturator colors:

#### 7) CONTRAINDICATION

No contraindication based on the patient population is known.

## 8) ADVERSE REACTIONS

None known.

## 9) WARNINGS AND PRECAUTIONS

For your safety, the safety of your patients and the safety of third parties, observe the safety notes below to minimize following risks to the Paste fillers instruments:

### User qualification and patient population

- The instructions for use are a component of the product and must be read carefully prior to use and be accessible at all times.
- The Paste fillers instruments may only be used in accordance with the intended use, any other type of use is not permitted.
- These instruments have not been tested on children, pregnant nor breastfeeding women.

### Infection / Toxic or allergic reaction hazards

- Paste fillers instruments must not be used on individuals with a known allergic sensitivity to stainless steel and silicone.
- For your own safety, wear personal protective equipment (masks, gloves, and goggles) during the treatment and reprocessing of the instrument.
- Do not use if non-sterile packaging is broken.
- Non-sterile instruments should be reprocessed before the first use and after each use.
- For reusable instruments, use approved cleaning and disinfecting agents:
  - Washing: Neodisher® Mediclean Forte (0.5 % concentration)
  - Thermal disinfection: Neodisher® Mediklar Special (0.03 % concentration)
 in accordance with the reprocessing instructions on FKG's website.
- Reprocess the instrument appropriately before disposal.

### Breakage / Inhalation or swallowing hazards

A damaged instrument or components could injure patients, users and third parties.

- Carefully read the label, marking on the packaging and this instruction for use to ensure the correct identification and use of device.
- Use a dental dam and a protective wire inside the safety hole when using manual endodontic instruments to avoid aspiration or ingestion by the patient.
- Take several radiographs from different angles to determine the anatomy of the root canals (length, width and curvature).
- Always inspect the instrument(s) before use, and discard (it) them if there are any visible defect(s).
- Cracks, deformations, signs of corrosion, loss in color or marking are signs that the device is no longer able to achieve the required performance level and should be discarded.
- Regularly inspect the instrument during use and discard if it shows any signs of distortion or wear.
- Do not exceed a sterilization temperature of 135 °C.

### Device packaging

- The packaging of the instruments should be carefully examined before opening (packaging integrity, no humidity, no visual degradation of instruments) to ensure that the packaging's integrity has not been compromised during transport and storage.

### Warnings and precautions for the (re)processing of instruments

- Paste fillers instruments that are provided non-sterile, must be reprocessed before the first use and between each use, following reprocessing instructions available on FKG's website.

Furthermore, it is the responsibility of the practitioner to always check its instruments before each use to identify possible signs of defects.

### In case of accident

- All serious events occurring in relation to the product must be reported to the manufacturer and the competent authority according to local regulations.

## 10) PROTOCOL FOR USE

- 1) Select a paste filler that is shorter than the root canal instrumentation length.
- 2) Apply the intracanal medication (e.g.,  $\text{Ca}(\text{OH})_2$ ) or obturation sealer (e.g., TotalFill) partially inside the root canal as recommended by the manufacturer.
- 3) Cover the paste filler tip with the medication or sealer.
- 4) Insert the paste filler into the root canal.
- 5) Operate the paste filler at a speed of 800 rpm and a torque limit of 1.5 Ncm (no torque limit is needed for Sensipast due to its automatic safety clutch).
- 6) Move the paste filler vertically up and down by about 4-6 mm, repeating the movement about 5 times per canal.
- 7) Remove the paste filler from the root canal while rotating.
- 8) Repeat the full protocol for each canal requiring treatment.

## 11) CONDITIONS AND TERM OF STORAGE, SHELF-LIFE

No specific shelf-life is indicated for the non-sterile instruments.

## 12) REPROCESSING PROTOCOL

Reprocessing must be applied to Paste fillers instruments before the first use and before each reuse. Instruments must be reprocessed following reprocessing instructions available on FKG's website.

Due to the design of the instruments and/or materials used and in the absence of contrary indications in the labeling or instructions for use of the device, the total number of uses is maximum 10 canals.

## 13) DISPOSAL

When an instrument reaches the end of its life, please dispose of it in accordance with the applicable laws and regulations for the safe disposal of sharp and contaminated instruments.

## 14) SYMBOLS

For explanation of symbols for IFUs and labels, please consult the document "Glossary of the symbols used by FKG Dentaire" available on FKG's website.

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