Multilingual Instruction for use XP-4D® Bio GP Cones & Absorbent Paper Points

Contents

EN	INSTRUCTION FOR USE.....	2
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА	7
CS	NÁVOD K POUŽITÍ.....	13
DA	BRUGSANVISNING	19
DE	GEBRAUCHSANWEISUNG	25
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	31
ES	INSTRUCCIONES DE USO	37
ET	KASUTUSJUHEND	43
FI	KÄYTTÖOHJEET	48
FR	NOTICE D'UTILISATION.....	53
HR	UPUTE ZA UPORABU	58
HU	HASZNÁLATI UTASÍTÁS.....	64
IT	ISTRUZIONI PER L'USO.....	69
LT	NAUDOJIMO INSTRUKCIJA	75
LV	LIETOŠANAS INSTRUKCIJAS	81
NL	GEBRUIKSAANWIJZING.....	86
NO	BRUKSANVISNING	92
PL	INSTRUKCJA UŻYWANIA	97
PT	INSTRUÇÕES DE UTILIZAÇÃO	102
RO	INSTRUCȚIUNI DE UTILIZARE.....	107
SK	NÁVOD NA POUŽITIE.....	112
SL	NAVODILA ZA UPORABO	118
SV	BRUKSANVISNING	124
TR	KULLANIM TALIMATLARI.....	129

FOR DENTAL USE ONLY

INSTRUCTIONS FOR USE – XP-4D® Bio GP Cones & Absorbent Paper Points

1) REFERENCES

XP-4D® Bio GP Cones						
References	08.600.03.025.RE	08.600.03.035.RE	08.600.03.045.RE	08.600.03.S25.RE		
Designation	STANDARD	MEDIUM	LARGE	Assortment S-L		
Color Code	Blue	Black	Yellow	Blue	Black	Yellow
Sterile	NO					
Reusable	NO					
Length	28 mm					
Bio GP Cones visual						
Correspondence with XP-4D® files	XP-4D® STANDARD 	XP-4D® MEDIUM 	XP-4D® LARGE 	XP-4D® STANDARD 	XP-4D® MEDIUM 	XP-4D® LARGE

XP-4D® Absorbent Paper Points						
References	08.610.03.030.BC	08.610.03.040.BC	08.610.03.050.BC	08.610.03.S30.BC		
Designation	STANDARD	MEDIUM	LARGE	Assortment S-L		
Color Code	Blue	Black	Yellow	Blue	Black	Yellow
Sterile	YES					
Reusable	NO					
Length	28 mm					
Absorbent Paper Points visual						
Correspondence with XP-4D® files	XP-4D® STANDARD 	XP-4D® MEDIUM 	XP-4D® LARGE 	XP-4D® STANDARD 	XP-4D® MEDIUM 	XP-4D® LARGE

2) PARAMETERS FOR USE

N/A

3) INTENDED USE

XP-4D® Bio GP Cones are intended to be used to fill a prepared root canal prior to tooth restoration.

XP-4D® Absorbent Paper Points are intended to be used to remove moisture from the root canal or to introduce medicaments. This product is a gamma-ray sterilized disposable medical device.

4) INDICATIONS FOR USE

XP-4D® Bio GP Cones are intended to be used to fill the root canal to prevent infection and instability of the root canal. It can be used with sealer or cement and it shall be used under the control of the dentist according to the clinical indication.

XP-4D® Absorbent Paper Points are designed to remove moisture from the root canal or to introduce medicaments.

XP-4D® Bio GP Cones & Absorbent Paper Points are intended to be used with the last XP-4D® file used for canal shaping.

XP-4D® Bio GP Cones & Absorbent Paper Points are intended for use in medical or hospital facilities, by qualified health professionals.

5) KEY PERFORMANCES

XP-4D® Bio GP Cones & Absorbent Paper Points are compliant with all applicable testing requirements of the ISO 13485 series of standards.

XP-4D® Bio GP Cones & Absorbent Paper size designation and length: Please contact the manufacturer for the size and length.

XP-4D® Bio GP Cones are radiopaque: Equivalent to 6mm aluminium stepwedge.

The SSCP (Summary of Safety and Clinical Performance) can be downloaded at <http://www.suredent.com>.

6) COMPOSITION

XP-4D® Bio GP Cones : Gutta Percha, Zinc Oxide, Zirconium Oxide, Barium Sulfate, Coloring Agent.

XP-4D® Absorbent Paper Points: Paper 100%

7) CONTRAINDICATION

Do not use XP-4D® Bio GP Cones for patient who is allergic to this product. Regarding XP-4D® Absorbent Paper Points, no contraindication based on the patient population is known.

8) ADVERSE REACTIONS

None known.

9) WARNINGS AND PRECAUTIONS

For your safety, the safety of your patients and the safety of third parties, observe the safety notes below to minimize the following risks to the XP-4D® Bio GP Cones & Absorbent Paper Points:

User qualification and patient population

- The instructions for use are a component of the product and must be read carefully prior to use and be accessible at all times.
- The XP-4D® Bio GP Cones & Absorbent Paper Points may only be used in accordance with the intended use, any other type of use is not permitted.
- Patient Population: All age group.

Infection / Toxic or allergic reaction hazards

- For your own safety, wear personal protective equipment (masks, gloves, and goggles) during the treatment.
- Do not use if sterile packaging seal is broken, damaged or wet.
- Do not use XP-4D® Bio GP Cones & Absorbent Paper Points after expiration date.
- Single use bio GP cones or absorbent paper points must not be retreated for subsequent use. This could compromise their integrity, and re-use of used absorbent paper point may cause inflammation.
- Use disinfected tools when bio GP cones or absorbent paper points are manipulated.

Breakage / Deterioration / Inhalation or swallowing hazards

A damaged bio GP cones or absorbent paper points or components could injure patients, users and third parties.

- Carefully read the label, marking on the packaging and this instruction for use to ensure the correct identification and use of device.
- Use a dental dam and a protective wire inside the safety hole when using bio GP cones or absorbent paper points to avoid aspiration or ingestion by the patient.
- Always inspect the bio GP cones absorbent paper point(s) before use, and discard (it) them if there are any visible defect(s).

Device packaging

- If the packaging has been opened, damaged, become wet or the expiry date has passed, the sterile state of the absorbent paper points inside the packaging is not guaranteed. Discard the device(s) that are single use.
- The packaging of the bio GP cones or absorbent paper points should be carefully examined before opening (packaging integrity, no humidity, no visual degradation of bio GP cones or absorbent paper points) to ensure that the packaging's integrity has not been compromised during transport and storage.

Furthermore, it is the responsibility of the practitioner to always check its devices before each use to identify possible signs of defects.

In case of accident

- All serious events occurring in relation to the product must be reported to the manufacturer and the competent authority according to local regulations.

10) PROTOCOL FOR USE

XP-4D® Bio GP Cones:

After shaping the root canal, use the proper size bio GP cones to fill.

XP-4D® Absorbent Paper Points:

After shaping the root canal, select the size that fits the file and dry the root canal, and use it with medication after.

11) CONDITIONS AND TERM OF STORAGE, SHELF-LIFE

The shelf-life is 4 years for XP-4D® Bio GP Cones & Absorbent Paper Points. The use-by date is indicated on the packaging. Sterility is preserved and guaranteed for the whole shelf life unless the packaging is broken, or the storage conditions are not respected. Store XP-4D® Bio GP Cones at room temperature away from direct sunlight, and store XP-4D® Absorbent Paper Points at room temperature.

12) REPROCESSING PROTOCOL

N/A

13) DISPOSAL

When a device reaches the end of its life, make sure that it is discarded in accordance with the applicable laws and regulations.

14) SYMBOLS

Information concerning the symbols used for labeling, packaging, and instructions for use:

Symbol		Designation	Description	Reference/ISO registration number
Identification of economic operators		Manufacturer	Identifies the manufacturer of the medical device (name and address). Note: If a date is next to this symbol, this corresponds to the date of manufacture. Date format: YYYY-MM-DD.	ISO 15223-1
		Authorized representative (authorized representative in the European Union)	Identifies the manufacturer's European authorized representative (name and address).	ISO 15223-1
		Authorized representative (authorized representative in Switzerland)	Identifies the manufacturer's Switzerland authorized representative (name and address).	ISO 15223-1
		Distributor	Identifies the distributor of the medical device (name and address).	ISO 15223-1
General symbols		Medical device	Indicates that this is a medical device.	(EU) 2017/745 Annex I, Article 23.2(q)
	 fkg.ch/ifu	Operating instructions	Indicates the need for the user to consult the operating instructions and/or the processing instructions for FKG devices, as well as this document, made available on FKG Dentaire's website.	ISO 15223-1
Product identification and other information specific to each batch		Catalog number	Indicates the manufacturer's catalog number so that the medical device can be identified.	ISO 15223-1
		Batch code	Indicates the manufacturer's batch code so that the batch can be identified.	ISO 15223-1


Symbol		Designation	Description	Reference/ISO registration number
		UDI code (Unique Device Identification)	Indicates the UDI code and UDI carrier specific to the device.	(EU) 2017/745 Annex I, Article 23.2(h)
		Quantity	Indicates the quantity contained in the packaging.	Not applicable
		Use-by date	Indicates the date after which the medical device is not to be used. Date format: YYYY-MM-DD.	ISO 15223-1
Other symbols relating to product identification and specifications		Do not reuse	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1
		Keep away from sunlight	Indicates a medical device that needs protection from light sources	ISO 15223-1
Information relating to sterility and (re)sterilization		Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation.	ISO 15223-1
		Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.	ISO 15223-1
		Non-sterile	Indicates a medical device that has not been subjected to a sterilization process.	ISO 15223-1
Other symbols		CE marking	European conformity marking for Class > I medical devices, accompanied by the identification number of the notified body. Notified body no. 0197: TÜV Rheinland LGA.	(EU) 2017/745 Annex V and 93/42/EC Annex XII

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Sure Dent Corporation
 #809,52, Sagimakgol-ro, Jungwon-gu,
 Seongnam-si, Gyeonggi-do, Korea
 E-mail : suredent@suredent.com

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OBELIS S.A
Obelis Swiss GmbH