## INSTRUCTIONS FOR USE – DirectFit™ Fiber Post

EN	Instructions for use .....	2
BG	Инструкции за употреба.....	8
CS	Návod k použití .....	14
DA	Brugsanvisning .....	20
DE	Gebrauchsanweisung .....	26
EL	Οδηγίες χρήσης.....	32
ES	Instrucciones de uso.....	38
ET	Kastutusjuhend .....	44
FI	Käyttöohjeet .....	49
FR	Notice d'utilisation .....	54
HR	Upute za uporabu.....	60
HU	Használati utasítás.....	66
IT	Istruzioni per l'uso.....	72
LT	Naudojimo instrukcija .....	77
LV	Lietošanas instrukcijas.....	82
NL	Gebruiksaanwijzing .....	87
NO	Bruksanvisning .....	92
PL	Instrukcja używania .....	97
PT	Instruções de utilização .....	103
RO	Instrucțiuni de utilizare .....	109
SK	Návod na použitie .....	114
SL	Navodila za uporabo .....	119
SV	Bruksanvisning .....	124
TR	Kullanım Talimatları .....	129

FOR DENTAL USE ONLY

## INSTRUCTIONS FOR USE – DirectFit™ Fiber Post

### DEFINITION

Intraradicular glass fiber posts for reinforcement of restorations and support of prosthetic crowns. Glass fiber posts enable rapid reconstruction of endodontically treated teeth and act in filling wide canals. The mechanical properties of the posts, being similar to dentin, prevent fracture and failure by minimizing the propagation and intensity of stresses in the root of the tooth.

The DirectFit™ Fiber Post is made of Zirconia AR-glass fibers that result in increased strength and radiopacity. The fibers are contained in a special translucent epoxy resin material which allows the DirectFit™ Fiber Post to be used with today's light or dual curing composite materials. The new DirectFit™ Fiber Post is factory silanated and exhibits a chemical adhesion to core composite material and bonded resin cements. The DirectFit™ Fiber Post's unique resin matrix allows for improved wetting and consequently, a higher ratio of fibers, thereby increasing the mechanical properties of the post.

### INTENDED USE

DirectFit™ Fiber Post is an endodontic intraradicular post, mainly composed of glass fiber and intended to be used in patients who require support for reinforcement of restorations and prosthetic crowns on endodontically treated teeth

### INDICATIONS

Glass fiber posts are indicated for retaining direct or indirect restorations in teeth with endodontic treatment.

### INTENDED USER

The device is intended to be used by trained dentists in endodontic and rehabilitative procedures.

### INTENDED ENVIRONMENT

Dental offices

### INTENDED POPULATION

DirectFit™ Fiber Post indicated for adult patients who require support and retention for prosthetic crowns and conventional restorations in cases of insufficient dental structure.

### CONTRAINDICATION

- Do not use DirectFit™ Fiber Posts in patients with reported sensitivity or allergy to any of its components.
- The use of DirectFit™ Fiber Posts is not indicated for extensive oral rehabilitations, patients with parafunctional habits, and teeth with less than 2 mm of remaining coronal structure.

- Do not use after the expiration date.
- Do not use the product if its packaging is damaged.

## COMPOSITION

ITEM	COMPOSITION
POST	Zirconia AR-glass fibers - 80% Epoxy Resin - 20%
STOP	Polyurethane

## WARNINGS

- Use protective goggles and a mask to prevent inhalation of particles generated during cutting of glass fiber posts.

*Note: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority health in which the user and/or patient is established.*

## PRECAUTIONS

- Use clinical tweezers to hold the DirectFit™ Fiber Post. Do not use hemostatic forceps or pliers, as they may crush the post.
- DirectFit™ Fiber Posts should only be cut with a diamond-tipped drill at high speed and under cooling. Do not use scissors or pliers. Use protective goggles and a mask to prevent inhalation of particles generated during cutting.
- According to endodontic care standards, always use rubber dam isolation when performing endodontic treatment or installing dental posts.
- DirectFit™ Fiber Posts should be fully covered by the cementing agent or composite resin. Exposed fibers absorb water, which alters the mechanical properties of the post.
- Use only sterile materials.
- Safety and efficacy of use have not been established in vulnerable populations.

## TECHNIQUE OF USE

### A. Root Canal Preparation

Fiber post installation should be done in rubber dam conditions

#### A.1. Preparation before the obturation step is performed

- Complete the root canal preparation.
- Select the post that fits the instrumented canal or select the post according to the canal diameter and depth using a periapical radiograph of the tooth. Keep in mind the amount of gutta percha you wish to be left in the canal (it's recommended to leave between 3 to 5 [mm] from apical foramen).
- Place the selected post in the root canal prior to obturation to confirm the fit and length of the post. Depending upon the desired length of the post in the canal, you can choose to go to a larger or smaller size post. It's better to have a little clearance than to have the post pinched in the root canal.

- To establish the length of the post, set the post cursor using the same reference point as used for the endodontic procedure. For a good hold, the length of the cut post should be at least half the length of the root canal and the emerging part should be at least half the height of the crown.
- Proceed to obturate the root canal system.

## A.2. Preparation after the obturation step is concluded

- Remove the required amount of gutta percha to accommodate the post, using a heated plugger or endodontic files (compatible with DirectFit™ Fiber Posts). Remove the gutta percha 1-2mm at a time, re-heating the plugger after each pass. It should take 3-5 passes depending on the length of the canal. As an alternative a powered heat source may be used. It's recommended to leave between 3 to 5 [mm] of gutta percha. Excessive removal of gutta percha may lead to failure.

DirectFit™ Fiber Posts			Compatible endodontic file	
Tip	Taper	Identification	Tip	Taper
050	0.04	Yellow	≤35	0.04
070		Red	40 - 45	
090		Blue	50 - 60	
110		Green	70 - 80	
050	0.06	Yellow	≤35	0.06
070		Red	40	
090		Blue	45 – 55	
100		Green	≥60	

- Prior to cementation, remove any obturation debris (gutta percha tags / sealer) which remain on the canal walls. This will ensure optimal bonding to the post and canal wall. The canal shape should not be altered.

## B. Preparation of DirectFit™ Fiber Post for cementation

- Clean the post and position it in the canal and take a radiography to check its adaptation and complete insertion.
- Slide the post cursor to the cutting area.
- Cut the post in the cervical portion with bur or disk with any kind of cooling to avoid resin degradation (diamond-tipped drills at high speed under water irrigation are recommended).
- Clean the post with alcohol and dry with air jets.
- Apply the adhesive system of choice to the post according to the manufacturer instructions.
- Do not polymerize the adhesive.

## C. Cementation

### C.1. Conventional dual-resin cement

- Rinse the canal with water jets and dry with air and absorbent paper points.
- Condition the canal with 37% phosphoric acid for 15 seconds.

- Rinse the applied acid for 30 seconds and dry the canal, leaving it slightly moist.

*Note that excess moisture can compromise the bonding strength between the canal wall and the resin cement.*

- Apply the adhesive system to the root canal according to the manufacturer's instructions.
- Remove excess adhesive with absorbent paper points.
- Do not polymerize the adhesive.
- Apply the conventional dual-resin cement to the post and fill the canal according to the cement instructions.
- Position the post in the canal, remove excesses, and perform photoactivation.

### **C.2. Self-adhesive dual-resin cement**

- Rinse the canal with water jets and dry with air and absorbent paper points.

*\*Due to the self-adhesive characteristics of the cement, the acid conditioning step should not be performed.*

- Apply the self-adhesive dual-resin cement to the post and fill the canal according to the cement instructions.
- Position the post in the canal, remove excesses, and perform photoactivation.

### **D. Fabrication of the coronal part of the core**

- Build the coronal part of the core with composite resin. The post should not be exposed outside the crown because of the risk of contamination.

### **STORAGE**

The product should be stored in a cool, dry place, protected from dirt.

### **DISPOSAL METHOD**

Dispose of the product according to national laws and recommendations from local authorities for proper disposal.

## SYMBOLS

For explanation of symbols for IFUs and labels, please consult the document “Glossary of the symbols used by FKG Dentaire” available on FKG’s website.

Symbols	Designation	Description	Reference/ISO registration number
	Manufacturer	Identifies the manufacturer of the medical device (name and address). Note: If a date is next to this symbol, this corresponds to the date of manufacture. Date format: YYYY-MM-DD.	ISO 15223-1 (ISO no. 7000-3082)
	CE marking	European conformity marking for Class I medical devices.	(EU) 2017/745 Annex V and 93/42/EEC Annex XII
 fkg.ch/ifu	Operating instructions	Indicates the need for the user to consult the operating instructions and/or the processing instructions for FKG devices, as well as this document, made available on FKG Dentaire's website.	ISO 15223-1 (ISO no. 7000-1641)
	Non-Sterile	Indicates a medical device that has not been subjected to a sterilization process.	ISO 15223-1 (ISO no. 7000-2609)
	Keep Dry	Indicates a medical device that needs to be protected from moisture	ISO 15223-1 (ISO no. 7000-0626)
	Keep away from sunlight	Indicates a medical device that needs protection from light sources	ISO 15223-1 (ISO no. 7000-0624)
	Do not reuse	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1 (ISO no. 7000-1051)
	Medical device	Indicates that this is a medical device.	(EU) 2017/745 Annex I, Article 23.2(q)
	Authorized representative (authorized representative in the European Union)	Identifies the manufacturer's European authorized representative (name and address).	ISO 15223-1
	Importer	Identifies the importer of the medical device (name and address).	ISO no. 7000-3725
	Catalog number	Indicates the manufacturer's catalog number so that the medical device can be identified.	ISO 15223-1 (ISO no. 7000-2493)

Symbols	Designation	Description	Reference/ISO registration number
	Batch code	Indicates the manufacturer's batch code so that the batch can be identified.	ISO 15223-1 (ISO no. 7000-2492)
	UDI code (Unique Device Identification)	Indicates the UDI code and UDI carrier specific to the device.	(EU) 2017/745 Annex I, Article 23.2(h)
	Quantity	Indicates the quantity contained in the packaging.	Not applicable
	Use-by date	Indicates the date after which the medical device is not to be used. Date format: YYYY-MM-DD.	ISO 15223-1 (ISO no. 7000-2607)
	Date of manufacture	Indicates the date on which the medical device is manufactured. Date format: YYYY-MM-DD.	ISO 15223-1 (ISO no. 7000-2497)

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