INSTRUCTIONS FOR USE – R-MOTION®

EN	Instructions for use	2
BG	Инструкции за употреба.....	7
CS	Návod k použití	12
DA	Brugsanvisning	17
DE	Gebrauchsanweisung	22
EL	Οδηγίες χρήσης	27
ES	Instrucciones de uso.....	32
ET	Kasutusjuhend.....	37
FI	Käyttöohjeet.....	42
FR	Notice d'utilisation	47
HR	Upute za uporabu	52
HU	Használati utasítás.....	57
IT	Istruzioni per l'uso.....	62
LT	Naudojimo instrukcija	67
LV	Lietošanas instrukcijas.....	72
NL	Gebruiksaanwijzing	77
NO	Bruksanvisning	82
PL	Instrukcja używania	87
PT	Instruções de utilização.....	92
RO	Instrucţiuni de utilizare	97
SK	Návod na použitie	102
SL	Navodila za uporabo	107
SV	Bruksanvisning	112
TR	Kullanım talimatları	117

FOR DENTAL USE ONLY

INSTRUCTIONS FOR USE – R-MOTION®

1) DESCRIPTION

R-MOTION®						
Glider C	Glider	20	25	30	40	50
15/.03v	15/.03v	20/.04	25/.06	30/.04	40/.04	50/.04
						
Non heat-treated	Heat-treated					

2) REFERENCES

References	Designation	Sizes	Sterile	Reusable	Length in mm		
					21	25	31
S1.7B0.00.0xx.FK	R-MOTION Glider C	015/.03v	YES	NO	X	X	X
	R-MOTION Glider	015/.03v	YES	NO	X	X	X
	R-MOTION 20	020/.04	YES	NO	X	X	X
	R-MOTION 25	025/.06	YES	NO	X	X	X
	R-MOTION 30	030/.04	YES	NO	X	X	X
	R-MOTION 40	040/.04	YES	NO	X	X	X
	R-MOTION 50	050/.04	YES	NO	X	X	X
99.700.80.S02.FK	Promo Set R-Motion	015/.03v – 025/.06 – 03/.04	YES	NO		X	

3) PARAMETERS FOR USE

R-Motion® instruments must only be used in **reciprocating** (counterclockwise cutting action) mode.
 R-Motion® instruments are compatible with endodontic motors compliant with ISO 1797.

4) INTENDED USE

R-Motion® endodontic instruments are intended for root canal treatment.

5) INDICATIONS FOR USE

R-Motion® endodontic instruments are indicated in case of pulp or periapical diseases that require root canal treatment.

R-Motion® Glider and R-Motion® Glider C instruments are designed for establishing a glide path of root canals in reciprocating mode.

R-Motion® 20/25/30/40/50 instruments are designed for shaping the root canals in reciprocating mode.


R-Motion® endodontic instruments are intended for use in medical or hospital facilities, by qualified health professionals.

6) KEY PERFORMANCES

R-Motion® instruments are compliant with all applicable testing requirements of the ISO 3630 series of standards.

7) COMPOSITION

The working part of these instruments is made of Nickel-Titanium and carries a shank (Crominox-plated maillechort). R-Motion® endodontic instruments carry a silicone endo stop to determine the working length of root canal preparation.

Example on 25 / 31 mm instruments

8) CONTRAINDICATION

No contraindication based on the patient population is known.

9) ADVERSE REACTIONS

None known.

10) WARNINGS AND PRECAUTIONS

For your safety, the safety of your patients and the safety of third parties, observe the safety notes below to minimize the following risks to R-Motion® endodontic instruments:

User qualification and patient population

- The instructions for use are a component of the product and must be read carefully prior to use and be accessible at all times.
- R-Motion® endodontic instruments may only be used in accordance with the intended use, any other type of use is not permitted.
- These instruments have not been tested on children, pregnant nor breastfeeding women.

Infection / Toxic or allergic reaction hazards

- R-Motion® endodontic instruments must not be used on individuals with a known allergic sensitivity to Nickel-Titanium and silicone.
- For your own safety, wear personal protective equipment (masks, gloves, and goggles) during the treatment and reprocessing of the instrument.
- Do not use if sterile packaging seal is broken, damaged or wet.
- Do not use the instruments after expiration date.

- Single use instruments must not be retreated for subsequent use. This could comprise their integrity.
- Reprocess the instrument appropriately before disposal.

Breakage / Deterioration / Inhalation or swallowing hazards

A damaged instrument or components could injure patients, users and third parties.

- Carefully read the label, marking on the packaging and this instruction for use to ensure the correct identification and use of device.
- Use a dental dam when using R-Motion® endodontic instruments to avoid inhalation or ingestion by the patient.
- Take several radiographs from different angles to determine the anatomy of the root canals (length, width and curvature).
- Always inspect the instrument(s) before use, and discard (it) them if there are any visible defect(s).
- Cracks, deformations, signs of corrosion, loss in color or marking are signs that the device is no longer able to achieve the required performance level and should be discarded.
- Regularly clean flutes and inspect instrument during use, and discard if it shows any signs of distortion or wear.
- If R-Motion® endodontic instrument is not progressing easily, withdraw it from the canal, clean it, check its cutting edges, then irrigate the canal and repeat the process.

Device packaging

- The packaging of the instruments should be carefully examined before opening (packaging integrity, no humidity, no visual degradation of instruments) to ensure that the packaging's integrity has not been compromised during transport and storage.
- If the packaging has been opened, damaged, become wet or the expiry date has passed, the sterile state of the instruments inside the packaging is not guaranteed. Discard the device(s).

Furthermore, it is the responsibility of the practitioner to always check its instruments before each use to identify possible signs of defects.

In case of accident

- All serious events occurring in relation to the product must be reported to the manufacturer and the competent authority according to local regulations.

11) PROTOCOL FOR USE


Figure 1

Figure 2

- 1) Create straight-line coronal and radicular access.
- 2) If necessary, use Barbed Broach to remove pulp tissue. *Barbed Broach must always enter the root canal gently.*
- 3) Use an ISO 010 hand file to explore the canal.
- 4) Determine the working length using a radiograph. In addition, electronical working length determination can be achieved.
- 5) Use the mechanized **R-Motion® Glider** instrument, applying gentle strokes of 2-3 mm to establish a glide path until the working length is reached. After 3 strokes (Fig. 1), remove and clean the instrument, then irrigate; recapitulate with an ISO 010 hand file. *Use of the **R-Motion® Glider C** is recommended for canals which require an instrument with more rigidity, and/or if the R-Motion® Glider is no longer progressing and can not reach the working length.*
- 6) Select the R-Motion® single-file shaping instrument which is best suited to the anatomy of the canal being treated. In most cases, and depending on your approach to apical preparation, the **R-Motion® 25** or **R-Motion® 30** will be appropriate. *If the radiographs do not allow you to clearly identify the anatomy of the canal, it can be gauged using hand files.*
- 7) Start the shaping procedure (Fig. 2), applying your usual irrigation protocol.
- 8) Use gentle inward-outward movements of 2-3 mm, applying very light apical pressure and allowing the instrument to progress passively along the canal. After 3 strokes, remove and clean the instrument, then irrigate; recapitulate with an ISO 010 hand file.
- 9) Repeat the previous step until the instrument has reached its working length. *If necessary, a brushing motion may be applied to remove coronal interferences and/or ensure uniform shaping of irregular canals. If the instrument is no longer progressing, make one or more passes with the R-Motion® 20 instrument, then resume using the R-Motion® 25 or the R-Motion® 30. Once the working length has been reached, remove the instrument to prevent over-enlargement of the foramen.*

Once all the canals have been shaped, proceed with the deep cleaning of the canals. *For an optimal result, the use of Finisher is recommended.*

Once the cleaning of the canals is completed, proceed to the next phase of treatment:

- 10) Obturate with gutta-percha and sealer. *TotalFill® BC Points™ and TotalFill® BC Sealer™ are recommended.*

12) CONDITIONS AND TERM OF STORAGE, SHELF-LIFE

The shelf-life is 5 years for sterile instruments. The use-by date is indicated on the packaging. Sterility is preserved and guaranteed for the whole shelf life unless the packaging is broken, or the storage conditions are not respected.

13) REPROCESSING PROTOCOL

R-Motion® instruments are for single use only. No reprocessing is allowed.

14) DISPOSAL

When a device reaches the end of its life, make sure that it is discarded in accordance with the applicable laws and regulations.

15) SYMBOLS

For explanation of symbols for IFUs and labels, please consult the document “Glossary of the symbols used by FKG Dentaire” available on FKG’s website.

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