INSTRUCTIONS FOR USE – RACE® EVO

EN	Instructions for use	2
BG	Инструкции за употреба.....	8
CS	Návod k použití	14
DA	Brugsanvisning	20
DE	Gebrauchsanweisung	26
EL	Οδηγίες χρήσης	32
ES	Instrucciones de uso.....	39
ET	Kasutusjuhend.....	45
FI	Käyttöohjeet.....	51
FR	Notice d'utilisation	57
HR	Upute za uporabu	63
HU	Használati utasítás.....	69
IT	Istruzioni per l'uso.....	75
LT	Naudojimo instrukcija	81
LV	Lietošanas instrukcijas.....	87
NL	Gebruiksaanwijzing	93
NO	Bruksanvisning	99
PL	Instrukcja użytkowania	105
PT	Instruções de utilização.....	111
RO	Instrucţiuni de utilizare	117
SK	Návod na použitie	123
SL	Navodila za uporabo	129
SV	Bruksanvisning	135
TR	Kullanım talimatları	141

FOR DENTAL USE ONLY

INSTRUCTIONS FOR USE – RACE® EVO

1) DESCRIPTION

Refill

RACE® EVO										
Opener	Glide Path									
25/.08	10/.02	15/.02	20/.02	10/.04	15/.04	15/.06				
										
Shaping										
30/.02	20/.04	25/.04	30/.04	35/.04	40/.04	50/.04	25/.06	30/.06	35/.06	40/.06
										
Recommended speed: 1'000 rpm										
Torque limit: 1.5 Ncm										

Sequences

RACE® EVO					
RACE EVO 4% sequence			RACE EVO 6% Sequence		
RE1	RE2	RE3	RE1	RE2	RE3
15/.04	25/.04	30/.04	15/.04	25/.04	25/.06
					
Recommended speed: 1'000 rpm					
Torque limit: 1.5 Ncm					

2) REFERENCES

References	Designation	Sizes	Sterile	Reusable*	Length in mm			
					19	21	25	31
S1.7B0.00.0xx.FK	RACE EVO 010/.02	010/.02	YES	YES		X	X	X
	RACE EVO 015/.02	015/.02	YES	YES		X	X	X
	RACE EVO 020/.02	020/.02	YES	YES		X	X	X
	RACE EVO 030/.02	030/.02	YES	YES		X	X	X
	RACE EVO 010/.04	010/.04	YES	YES		X	X	X
	RACE EVO 015/.04	015/.04	YES	YES		X	X	X
	RACE EVO 020/.04	020/.04	YES	YES		X	X	X
	RACE EVO 025/.04	025/.04	YES	YES		X	X	X
	RACE EVO 030/.04	030/.04	YES	YES		X	X	X
	RACE EVO 035/.04	035/.04	YES	YES		X	X	X
	RACE EVO 040/.04	040/.04	YES	YES		X	X	X
	RACE EVO 050/.04	050/.04	YES	YES		X	X	X
	RACE EVO 015/.06	015/.06	YES	YES		X	X	X
	RACE EVO 020/.06	020/.06	YES	YES		X	X	X
	RACE EVO 025/.06	025/.06	YES	YES		X	X	X
	RACE EVO 030/.06	030/.06	YES	YES		X	X	X
	RACE EVO 035/.06	035/.06	YES	YES		X	X	X
	RACE EVO 040/.06	040/.06	YES	YES		X	X	X
	RACE EVO 025/.08	025/.08	YES	YES	X			
99.700.80.S03.FK	Promo Set RACE EVO	015/.04 - 025/.04 - 030/.04 - 025/.06	YES	YES			X	

* These devices can be used for a maximum of 8 canals and can be reprocessed up to 10 times.

3) PARAMETERS FOR USE

Recommended speed: 1'000 rpm continuous rotation

Torque limit: 1.5 Ncm

RACE® EVO instruments are compatible with endodontic motors compliant with ISO 1797.

4) INTENDED USE

RACE® EVO endodontic instruments are intended for root canal treatment.

5) INDICATIONS FOR USE

RACE® EVO endodontic instruments are indicated in case of pulp or periapical diseases that require root canal treatment.

All instruments are designed for preparation and shaping of root canals.


RACE® EVO endodontic instruments are intended for use in medical or hospital facilities, by qualified health professionals.

6) KEY PERFORMANCES

RACE® EVO instruments are compliant with all applicable testing requirements of the ISO 3630 series of standards.

7) COMPOSITION

The working part of these instruments is made of Nickel-Titanium and carries a shank (Crominox-plated maillechort). RACE® EVO endodontic instruments carry a silicone endo stop to determine the working length of root canal preparation.

Example on 25 / 31 mm instruments

Identification of Taper and Tip diameter using dual-color coding on the shank.

Identification of the working part length with silicone endo stop colors:

8) CONTRAINDICATION

No contraindication based on the patient population is known.

9) ADVERSE REACTIONS

None known.

10) WARNINGS AND PRECAUTIONS

For your safety, the safety of your patients and the safety of third parties, observe the safety notes below to minimize the following risks to RACE® EVO endodontic instruments:

User qualification and patient population

- The instructions for use are a component of the product and must be read carefully prior to use and be accessible at all times.
- RACE® EVO endodontic instruments may only be used in accordance with the intended use, any other type of use is not permitted.
- These instruments have not been tested on children, pregnant nor breastfeeding women.

Infection / Toxic or allergic reaction hazards

- RACE® EVO endodontic instruments must not be used on individuals with a known allergic sensitivity to Nickel-Titanium and silicone.
- For your own safety, wear personal protective equipment (masks, gloves, and goggles) during the treatment and reprocessing of the instrument.
- Do not use if sterile packaging seal is broken, damaged or wet.
- Do not use the instruments after expiration date.
- For reusable instruments, use approved cleaning and disinfecting agents:
 - Washing: Neodisher® Mediclean Forte (0.5 % concentration)
 - Thermal disinfection: Neodisher® Mediklar Special (0.03 % concentration)
 in accordance with the reprocessing instructions (see section 13).
- Reprocess the instrument appropriately before disposal.

Breakage / Deterioration / Inhalation or swallowing hazards

A damaged instrument or components could injure patients, users and third parties.

- Carefully read the label, marking on the packaging and this instruction for use to ensure the correct identification and use of device.
- Use a dental dam when using RACE® EVO endodontic instruments to avoid inhalation or ingestion by the patient.
- Take several radiographs from different angles to determine the anatomy of the root canals (length, width and curvature).
- Always inspect the instrument(s) before use, and discard (it) them if there are any visible defect(s).
- Cracks, deformations, signs of corrosion, loss in color or marking are signs that the device is no longer able to achieve the required performance level and should be discarded.
- Regularly clean flutes and inspect instrument during use, and discard if it shows any signs of distortion or wear.
- If RACE® EVO endodontic instrument is not progressing easily, withdraw it from the canal, clean it, check its cutting edges, then irrigate the canal and repeat the process.
- Do not exceed a sterilization temperature of 135°C.

Device packaging

- The packaging of the instruments should be carefully examined before opening (packaging integrity, no humidity, no visual degradation of instruments) to ensure that the packaging's integrity has not been compromised during transport and storage.
- If the packaging has been opened, damaged, become wet or the expiry date has passed, the sterile state of the instruments inside the packaging is not guaranteed. Discard the device(s) that are single use and consult the reprocessing instructions available on FKG's website for multiple use device(s).

Warnings and precautions for the (re)processing of instruments

- RACE® EVO endodontic instruments that are provided sterile, may be reprocessed for multi-usage following reprocessing instructions available on FKG's website

Furthermore, it is the responsibility of the practitioner to always check its instruments before each use to identify possible signs of defects.

In case of accident

- All serious events occurring in relation to the product must be reported to the manufacturer and the competent authority according to local regulations.

11) PROTOCOL FOR USE


Figure 1


Figure 2


Figure 3

- 1) Create straight-line coronal and radicular access.
- 2) If necessary, use Barbed Broach to remove pulp tissue. *Barbed Broach must always enter the root canal gently.*
- 3) Use an ISO 010 hand file to explore the canal.
- 4) Determine the working length using a radiograph. In addition, electronic working length determination can be achieved.
- 5) Use the RACE® EVO instruments in the order predefined by the chosen sequence. Start with glide path, followed by shaping of the canal.
- 6) RACE® EVO instruments are to be used with gentle strokes of 2-3 mm, applying very light apical pressure and allowing the instrument to progress passively along the canal. After 3 strokes, remove and clean the instrument, then irrigate ; recapitulate with an ISO 010 hand file.
- 7) Use the mechanized **RE1** instrument (Fig. 1) to perform glide path until the working length is reached.
- 8) Start the shaping procedure, applying your usual irrigation/disinfection protocol.
- 9) Use the **RE2** (Fig. 2) until the working length is reached.
- 10) Use the **RE3** (Fig. 3) until the instrument can be inserted to its working length. *If necessary, a brushing motion may be applied to remove coronal interferences and/or ensure uniform shaping of irregular canals.* Once the working length has been reached, remove the instrument to prevent over-enlargement of the foramen.
- 11) If the preparation is not sufficient for the anatomy of the canal being treated, then continue shaping the canal using the larger sizes of RACE® EVO instruments. *We recommend that **40/.04** or **50/.04** instruments are used for shaping larger canals, following the **RACE EVO 4%** sequence. We recommend that **30/.06** or **35/.06** instruments are used for shaping larger canals, following the **RACE EVO 6%** sequence.*

Once all the canals have been shaped, proceed with the deep cleaning of the canals. *For an optimal result, the use of Finisher is recommended.*

Once the cleaning of the canals is completed, proceed to the next phase of treatment:

12) Obturate with gutta-percha and sealer. *TotalFill® BC Points™ and TotalFill® BC Sealer™ are recommended.*

12) CONDITIONS AND TERM OF STORAGE, SHELF-LIFE

The shelf-life is 5 years for sterile instruments. The use-by date is indicated on the packaging.

Sterility is preserved and guaranteed for the whole shelf life unless the packaging is broken, or the storage conditions are not respected.

13) REPROCESSING PROTOCOL

Instruments must be reprocessed following reprocessing instructions available on FKG's website.

Due to the design of the instruments and/or materials used and in the absence of contrary indications in the labeling or instructions for use of the device, the number of uses is maximum 8 canals and can be reprocessed up to 10 times.

14) DISPOSAL

When a device reaches the end of its life, make sure that it is discarded in accordance with the applicable laws and regulations.

15) SYMBOLS

For explanation of symbols for IFUs and labels, please consult the document "Glossary of the symbols used by FKG Dentaire" available on FKG's website.

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2) REFERÊNCIAS

Referências	Designação	Tamanhos	Estéril	Reutilizável*	Comprimento em mm			
					19	21	25	31
S1.7B0.00.0xx.FK	RACE EVO 010/.02	010/.02	SIM	SIM		X	X	X
	RACE EVO 015/.02	015/.02	SIM	SIM		X	X	X
	RACE EVO 020/.02	020/.02	SIM	SIM		X	X	X
	RACE EVO 030/.02	030/.02	SIM	SIM		X	X	X
	RACE EVO 010/.04	010/.04	SIM	SIM		X	X	X
	RACE EVO 015/.04	015/.04	SIM	SIM		X	X	X
	RACE EVO 020/.04	020/.04	SIM	SIM		X	X	X
	RACE EVO 025/.04	025/.04	SIM	SIM		X	X	X
	RACE EVO 030/.04	030/.04	SIM	SIM		X	X	X
	RACE EVO 035/.04	035/.04	SIM	SIM		X	X	X
	RACE EVO 040/.04	040/.04	SIM	SIM		X	X	X
	RACE EVO 050/.04	050/.04	SIM	SIM		X	X	X
	RACE EVO 015/.06	015/.06	SIM	SIM		X	X	X
	RACE EVO 020/.06	020/.06	SIM	SIM		X	X	X
	RACE EVO 025/.06	025/.06	SIM	SIM		X	X	X
	RACE EVO 030/.06	030/.06	SIM	SIM		X	X	X
	RACE EVO 035/.06	035/.06	SIM	SIM		X	X	X
	RACE EVO 040/.06	040/.06	SIM	SIM		X	X	X
	RACE EVO 025/.08	025/.08	SIM	SIM	X			
99.700.80.S03.FK	Promo Set RACE EVO	015/.04 - 025/.04 - 030/.04 - 025/.06	SIM	SIM			X	

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