INSTRUCTIONS FOR USE – BC SEALER ION+

EN	Instructions for use	3
BG	Инструкции за употреба	10
CS	Návod k použití	18
DA	Brugsanvisning	26
DE	Gebrauchsanweisung	34
EL	Οδηγίες χρήσης.....	42
ES	Instrucciones de uso.....	50
ET	Kastutusjuhend	57
FI	Käyttöohjeet	64
FR	Notice d'utilisation	71
HR	Upute za uporabu.....	79
HU	Használati utasítás.....	86
IT	Istruzioni per l'uso.....	94
LT	Naudojimo instrukcija	102
LV	Lietošanas instrukcijas.....	110
NL	Gebruiksaanwijzing	118
NO	Bruksanvisning	126
PL	Instrukcja używania	133
PT	Instruções de utilização	141
RO	Instrucțiuni de utilizare	149

SK	Návod na použitie	157
SL	Navodila za uporabo	165
SV	Bruksanvisning	173
TR	Kullanım Talimatları	180

FOR DENTAL USE ONLY

INSTRUCTIONS FOR USE – BC SEALER ION+

INTENDED USE

BC SEALER ION+ is a bioceramic root canal sealer, intended to be used exclusively by dentists during endodontic obturation procedures. It is designed to fill and seal previously cleaned and shaped root canal systems, with the goal of preventing reinfection and ensuring long-term obturation. The device is provided in a ready-to-use presentation.

INDICATIONS

- Obturation of root canals of permanent teeth.
- Used during the treatment of internal resorption.

INTENDED POPULATION

Patients presenting permanent teeth with closed apices, diagnosed with pulpal and periapical disease.

INTENDED USER

BC SEALER ION+ is intended for use by dental surgeons qualified to perform endodontic procedures.

INTENDED ENVIRONMENT

Dental offices/clinics.

CONTRAINDICATIONS

- Do not use in patients with reported sensitivity to the components of the formula.
- Do not use the product after the expiration date.
- Do not use the product if the packaging is damaged.
- Do not reuse the applicator tips.

COMPOSITION

Calcium Magnesium Silicate, Calcium Sulfate, Potassium Sulfate, Calcium Oxide, Zirconium Oxide, Silicon Dioxide, and Dispersing Agent.

TECHNICAL INFORMATION

- Biocompatible product.
- Setting time: 30 – 240 minutes after insertion into the canal.

IMPORTANT: Dry with paper points, but without excessively drying the canal. Moisture from the dentinal tubules will initiate the setting reaction of the material.

- Ready to use, requires no mixing.
- High alkalinity.
- High radiopacity.
- Can be used with single cone and thermoplastic obturation techniques, in addition to conventional techniques.

WARNINGS

- Do not use BC SEALER ION+ for canal filling without gutta-percha cones.
- Avoid contact with eyes or skin. In case of contact, rinse immediately with water.
- Avoid contact with the oral mucosa. In case of contact, rinse with water and avoid ingestion of the product. If any sensitivity persists, seek medical attention immediately.
- The safety and efficacy of use have not been established in vulnerable populations.
- Keep out of reach of children.

NOTE: Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent health authority in which the user and/or patient is established.

PRECAUTIONS

- Use protective eyewear, mask, and gloves when handling the product. In case of contact with eyes or skin, wash with water.
- Use absolute isolation (rubber dam) for all indications.
- The product is moisture-sensitive; therefore, close the package properly with adequate pressure to prevent the product from drying out.
- Use the product only after controlling the infection, as the material's setting does not occur at this stage due to the acidic pH of the medium.
- Avoid causing excessive extrusion of the product into the periapical region.
- Disinfect the applicator tip prior to use. After use, the tip must be discarded.
- Clean and disinfect the syringe for subsequent uses.

- Return the properly closed syringe, along with the silica sachet, inside the packaging, which must be correctly closed and stored according to the "Storage" section.
- For subsequent use, check if the silica sachet inside the envelope remains blue in color.

USAGE TECHNIQUE USE OF APPLICATOR TIPS

1. It is recommended to use a disposable plastic cover on the syringe to avoid cross-contamination.
2. Confirm the flow of material from the syringe prior to direct application into the canal.
3. Properly position the tip and syringe assembly in the canal and press the plunger with light force to avoid excessive output of the product.
4. Before closing the syringe, retract the plunger to stop the outflow.
5. Clean the excess material from the syringe tip and close with adequate pressure after use to prevent the entry of moisture, which causes the product to dry out.
6. Clean and disinfect the syringe for subsequent uses to avoid cross-contamination. In case of contamination with saliva or blood, discard it.

A. OBTURATION OF ROOT CANALS OF PERMANENT TEETH TRADITIONAL TECHNIQUE - LATERAL COMPACTION:

1. Anesthetize, install absolute isolation, and then perform biomechanical preparation of the canal.
2. Remove the sealing cap and attach the applicator tip to the syringe.
3. Dry the canal only with paper points without causing excessive desiccation.
4. Apply BC SEALER ION+ by gradually retracting the syringe to obtain a complete filling of the canals.

Important: Press the plunger with light force to avoid excessive output of the product.

5. Introduce the master gutta-percha cone into the canal according to the last file used, followed by accessory cones, all coated with BC SEALER ION+.
6. Take a radiograph to verify the correct filling of the canal.
7. Cut the cone at the desired height with heated instruments followed by vertical compaction.
8. Remove excess material from the canal walls with water, perform coronal sealing and restoration.

Important:

- BC SEALER ION+ does not undergo setting contraction.
- BC SEALER ION+ can be used in conventional techniques, single cone, and thermoplastic obturation.
- Use conventional gutta-percha removal techniques.

SINGLE CONE TECHNIQUE

1. Anesthetize, install absolute isolation, and then perform biomechanical preparation of the canal.
2. Remove the sealing cap and attach the applicator tip to the syringe.
3. Dry the canal only with paper points without causing excessive desiccation.
4. Apply BC SEALER ION+ by gradually retracting the syringe to obtain a complete filling of the canals.
 - **IMPORTANT:** Press the plunger with light force to avoid excessive output of the product.
5. Introduce the gutta-percha cone.
6. Take a radiograph to verify the correct filling of the canal.
7. Cut the cone at the desired height with heated instruments followed by vertical compaction.
8. Remove excess material from the canal walls with water and appropriate irrigating solutions for this purpose, perform coronal sealing and restoration.

B. USE DURING INTERNAL RESORPTION TREATMENT

1. Anesthetize and install absolute isolation.
2. Remove the granulation tissue from the resorption area with sharp curettes.
3. Neutralize the environment with calcium hydroxide paste or another intracanal medicament of your preference.
4. Remove the calcium hydroxide or chosen intracanal medicament at the next session.
5. Dry the canal only with paper points without causing excessive desiccation.
6. Insert BC SEALER ION+ with the applicator tip throughout the canal, prioritizing the resorption site.
7. Introduce accessory gutta-percha cones according to the lateral compaction technique.
8. Cut the cones at the upper part of the resorption with heated instruments. Perform vertical compaction for better flow of the cement into the resorption site.
9. Take a radiograph to verify the correct filling of the resorption site and the canal.
10. Obturate the remainder of the canal with gutta-percha.
11. Perform coronal sealing with glass ionomer or another material of your preference and restore the tooth.
12. Radiograph and perform monitoring (follow-up) for at least two years.

REMOVAL OF OBTURATION

1. Use conventional gutta-percha removal techniques.

STORAGE

- Keep the product in a dry and ventilated place between 15 and 30°C (59 - 86°F) and relative air humidity below 60%.
- Do not store in a refrigerator.

DISPOSAL METHOD

Dispose of the product in accordance with national laws and local authority recommendations for proper disposal. Inadequate disposal may present biological or physical risks; handle and dispose of as hazardous biological material if contaminated.

SYMBOLS

For explanation of symbols for IFUs and labels, please consult the document “Glossary of the symbols used by FKG Dentaire” available on FKG’s website.

Symbols	Designation	Description	Reference/ISO registration number
	Manufacturer	Identifies the manufacturer of the medical device (name and address). Note: If a date is next to this symbol, this corresponds to the date of manufacture. Date format: YYYY-MM-DD.	ISO 15223-1 (ISO no. 7000-3082)
	CE marking	European conformity marking for Class > I medical devices, accompanied by the identification number of the notified body. Notified body no. 2460: DNV Product Assurance AS	(EU) 2017/745 Annex V and 93/42/EEC Annex XII
 fkg.ch/ifu	Operating instructions	Indicates the need for the user to consult the operating instructions and/or the processing instructions for FKG devices, as well as this document, made available on FKG Dentaire's website.	ISO 15223-1 (ISO no. 7000-1641)
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions to be taken that cannot be presented on the medical device itself.	ISO 9687 (ISO no. 7000-0434A)

Symbols	Designation	Description	Reference/ISO registration number
	Keep Dry	Indicates a medical device that needs to be protected from moisture	ISO 15223-1 (ISO no. 7000-0626)
	Temperature Limit	Indicates the temperature limits to which the medical device can be safely exposed	ISO 15223-1 (ISO no. 7000-0632)
	Do not reuse	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1 (ISO no. 7000-1051)
	Medical device	Indicates that this is a medical device.	(EU) 2017/745 Annex I, Article 23.2(q)
	Authorized representative (authorized representative in the European Union)	Identifies the manufacturer's European authorized representative (name and address).	ISO 15223-1
	Catalog number	Indicates the manufacturer's catalog number so that the medical device can be identified.	ISO 15223-1 (ISO no. 7000-2493)
	Batch code	Indicates the manufacturer's batch code so that the batch can be identified.	ISO 15223-1 (ISO no. 7000-2492)
	UDI code (Unique Device Identification)	Indicates the UDI code and UDI carrier specific to the device.	(EU) 2017/745 Annex I, Article 23.2(h)
	Quantity	Indicates the quantity contained in the packaging.	Not applicable
	Use-by date	Indicates the date after which the medical device is not to be used.	ISO 15223-1 (ISO no. 7000-2607)

Symbols	Designation	Description	Reference/ISO registration number
		Date format: YYYY-MM-DD.	
	Date of manufacture	Indicates the date on which the medical device is manufactured. Date format: YYYY-MM-DD.	ISO 15223-1 (ISO no. 7000-2497)
	Authorized representative (authorized representative in the Switzerland)	Identifies the manufacturer's Swiss authorized representative (name and address).	ISO 15223-1
	Non-sterile	Indicates that the device that is normally provided sterile in the same or similar packaging has not been sterilized.	ISO 15223-1 (ISO no. 7000-2609)
	Keep away from sunlight	Indicates a medical device that needs to be protected from sunlight	ISO 15223-1 (ISO no. 7000-0624)

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