

## Speed & torque

Product	Speed - minimum	Speed - optimal	Torque
Own sequence	400	600	1,5
PreRace	400	600	1,5
ScoutRace	600	800	1
Race ISO10	600	800	1,5
BT-Apisafe	600	800	1
BT-Race	800	1'000	1,5
BT-Race XL	600	800	1,5
iRace	400	600	1,5
iRace Plus	400	600	1,5
BioRace	600	600	1,5
BioRace Extended	600	600	1,5
EasyRace	400	600	1,5
XtremRace	400	600	1,5
D-Race DR1	800	1'000	1,5
D-Race DR2	400	600	1,5
XP-endo Finisher	800	800-1'000	1