


# Change the Game



Catalogue



For over 90 years, we have striven to deliver Swiss quality, precision and peace of mind to endo specialists. Our unique approach is driven by the quest to provide the best dental care and our conviction that the search for perfection is never over. Like you, we seek to advance the frontiers of the possible by evolving through research, refinement and innovation.  
Ready to change the game?

# Change the Game

We relentlessly challenge the status quo and offer constant innovation and high-quality solutions to preserve the tooth and improve Endodontics



Innovation

Engagement

Responsibility

## Table of contents

<b>A Equipment</b> .....	<b>4-5</b>
<b>01 Motor and apex locator</b> .....	<b>4</b>
Router X3000 .....	5
Router Universal .....	5
<b>B Files</b> .....	<b>6-17</b>
<b>01 Opening</b> .....	<b>6</b>
Gates & XS Gates .....	7
Peeso .....	7
RACE EVO .....	7
PreRace .....	7
<b>02 Glide path</b> .....	<b>8</b>
XP-endo Rise .....	9
R-Motion .....	9
RACE EVO .....	9
Race .....	9
<b>03 Shaping</b> .....	<b>10</b>
XP-endo Rise .....	11
R-Motion .....	11
RACE EVO .....	12
<b>04 Cleaning</b> .....	<b>14</b>
Finisher .....	15
<b>05 Retreatment</b> .....	<b>16</b>
XP-endo Rise .....	17
D-Race .....	17
<b>C Hand Files</b> .....	<b>18-23</b>
<b>01 Hand Files</b> .....	<b>18</b>
K-Files .....	19
K-Files Flex .....	20
H-Files .....	21
Reamers .....	22
F-Search .....	23
Barbed Broaches .....	23
<b>D Obturation</b> .....	<b>24-29</b>
<b>01 Bioceramic obturation</b> .....	<b>24</b>
TotalFill .....	25
BC Universal .....	26
<b>02 Paper Points &amp; Gutta Percha</b> .....	<b>27</b>
<b>03 Paste fillers</b> .....	<b>28</b>
Lentulo & Spring Paste Fillers .....	28
Sensipast .....	28
<b>04 Finger Spreaders &amp; Pluggers</b> .....	<b>29</b>
<b>E Accessories</b> .....	<b>30-31</b>
<b>F Restoration</b> .....	<b>32-33</b>
<b>G Legend</b> .....	<b>34</b>



# A. Equipment



# 01

## Motor and apex locator



### Rooter® X3000



Name	Ref.
Rooter X3000	08.971.00.000.FK
Rooter X3000 - Motor Handpiece	08.971.00.001.FK
Rooter X3000 - Charging Base	08.971.00.002.FK
Rooter X3000 - Contra-angle	08.971.00.003.FK
Rooter X3000 - Spray Nozzle	08.971.00.004.FK
Rooter X3000 - Measuring Wire - Mini USB	08.971.00.005.FK
Rooter X3000 - File Clip	08.971.00.006.FK
Rooter X3000 - Lip Hook	08.971.00.007.FK
Rooter X3000 - Touch Probe	08.971.00.008.FK
Rooter X3000 - Protective Cover	08.971.00.009.FK
Rooter X3000 - O-ring	08.971.00.010.FK
Rooter X3000 - AC-Adapter	08.971.00.011.FK
Rooter X3000 - Demo Version	08.971.00.012.FK
Rooter X3000 - Lithium-ion Battery	08.971.00.013.FK
Rooter X3000 - Measuring Wire - USB-C	08.971.00.014.FK



### Rooter® Universal




Name	Ref.
Rooter Universal	08.970.00.000.FK
Rooter Universal - Motor Handpiece	08.970.00.001.FK
Rooter Universal - Battery Charger	08.970.00.002.FK
Rooter Universal - AC-Adapter	08.970.00.003.FK
Rooter Universal - Tester	08.970.00.004.FK
Rooter Universal - Test Wire A	08.970.00.005.FK
Rooter Universal - Test Wire B	08.970.00.006.FK
Rooter Universal - File Clip	08.970.00.007.FK
Rooter Universal - Lip Hook	08.970.00.008.FK
Rooter Universal - Lighting Device	08.970.00.009.FK
Rooter Universal - Protective Sleeve	08.970.00.010.FK
Rooter Universal - Contra-angle	08.970.00.011.FK
Rooter Universal - Wrench (Contra-angle tool)	08.970.00.012.FK
Rooter Universal - Spray Nozzle	08.970.00.013.FK
Rooter Universal - AC-Adapter AU/NZ	08.970.00.014.FK
Rooter Universal - AC-Adapter UK	08.970.00.015.FK
Rooter Universal - Dust Cap	08.970.00.016.FK
Rooter Universal - Silicone Plug	08.970.00.017.FK
Rooter Universal - Demo Version	08.970.00.018.FK
Rooter Universal - Lithium-ion Battery	08.970.00.019.FK




# B. Files

## 01

### Opening




#### Gates & XS Gates



Length (mm)	ISO	Name	Plastic Box 6 pcs	Assortment 6 pcs
15	50	XS Gates 1	03.140.10.601.GC	XS Gates 1-6 03.140.10.S16.GC
	70	XS Gates 2	03.140.10.602.GC	
	90	XS Gates 3	03.140.10.603.GC	
	110	XS Gates 4	03.140.10.604.GC	
	130	XS Gates 5	03.140.10.605.GC	
	150	XS Gates 6	03.140.10.606.GC	
19	50	Gates 1	03.140.19.601.GC	Gates 1-6 03.140.19.S16.GC
	70	Gates 2	03.140.19.602.GC	
	90	Gates 3	03.140.19.603.GC	
	110	Gates 4	03.140.19.604.GC	
	130	Gates 5	03.140.19.605.GC	
	150	Gates 6	03.140.19.606.GC	




#### Peeso



Length (mm)	ISO	Name	Plastic Box 6 pcs	Assortment 6 pcs
19	70	Peeso 1	03.132.19.601.GC	Peeso 1-6 03.132.19.S16.GC
	90	Peeso 2	03.132.19.602.GC	
	110	Peeso 3	03.132.19.603.GC	
	130	Peeso 4	03.132.19.604.GC	
	150	Peeso 5	03.132.19.605.GC	
	170	Peeso 6	03.132.19.606.GC	





#### RACE® EVO



Length (mm)	ISO	Taper	Blister 6 pcs
19	25	0.08	S1.7B0.00.0GL.FK



#### PreRace




Length (mm)	ISO	Taper	Blister 6 pcs
19	30	0.06	S1.7B0.00.0DD.FK
	35	0.08	S1.7B0.00.0DA.FK
	40	0.10	S1.7B0.00.0DF.FK



# B. Files

## 02

### Glide path



#### XP-endo® Rise Glider



⌚ NITI STERILE R

Length (mm)	ISO	Taper	Blister 6 pcs
21	15 O	0.04	S1.7B0.00.0HN.FK
25			S1.7B0.00.0HO.FK
31			S1.7B0.00.0HP.FK



#### R-Motion®



⌚ NITI STERILE R

Length (mm)	ISO	Taper	Name	Blister 3 pcs
21	15 O	0.03v	R-Motion Glider	S1.7B0.00.0GN.FK
			R-Motion Glider C	S1.7B0.00.0HT.FK
25			R-Motion Glider	S1.7B0.00.0GN.FK
			R-Motion Glider C	S1.7B0.00.0HU.FK
31			R-Motion Glider	S1.7B0.00.0GO.FK
			R-Motion Glider C	S1.7B0.00.0HV.FK



#### RACE® EVO



⌚ NITI STERILE R

Length (mm)	ISO	Taper	Blister 6 pcs
21	10	0.02	S1.7B0.00.0GV.FK
	15 O		S1.7B0.00.0GX.FK
	20		S1.7B0.00.0GZ.FK
	10	0.04	S1.7B0.00.0HD.FK
	15 O		S1.7B0.00.0DG.FK
	15 O		S1.7B0.00.0HJ.FK
25	10	0.02	S1.7B0.00.0FU.FK
	15 O		S1.7B0.00.0FV.FK
	20		S1.7B0.00.0FW.FK
	10	0.04	S1.7B0.00.0FY.FK
	15 O		S1.7B0.00.0FZ.FK
	15 O		S1.7B0.00.0GG.FK
31	10	0.02	S1.7B0.00.0GW.FK
	15 O		S1.7B0.00.0GY.FK
	20		S1.7B0.00.0HA.FK
	10	0.04	S1.7B0.00.0HE.FK
	15 O		S1.7B0.00.0DH.FK
	15 O		S1.7B0.00.0HK.FK



#### Race



⌚ NITI STERILE R

Length (mm)	ISO	Taper	Blister 6 pcs	
21	10	0.02	S1.7B0.00.0AA.FK	
	15 O		S1.7B0.00.0AB.FK	
	20		S1.7B0.00.0AC.FK	
25	10	0.04	S1.7B0.00.0BH.FK	
	10		0.02	S1.7B0.00.0AL.FK
	15 O			S1.7B0.00.0AM.FK
31	20	0.04		S1.7B0.00.0AN.FK
	10		S1.7B0.00.0BP.FK	
	10		0.02	S1.7B0.00.0AW.FK
	15 O	S1.7B0.00.0AX.FK		
	20	S1.7B0.00.0AY.FK		
		10	0.04	S1.7B0.00.0BX.FK






# B. Files

## 03

### Shaping



#### XP-endo® Rise



NITI STERILE R

Length (mm)	ISO	Taper	Name	Blister 3 pcs	Blister 6 pcs	Sequences 6pcs
21	15 O	0.04	XP-endo Rise Glider	-	S1.7B0.00.0HN.FK	3 x (XP-RG + XP-RS) S1.XB0.00.SAI.FK
	30 ●		XP-endo Rise Shaper	S1.XB0.00.0AU.FK	S1.XB0.00.0AR.FK	
25	15 O	0.04	XP-endo Rise Glider	-	S1.7B0.00.0HO.FK	3 x (XP-RG + XP-RS) S1.XB0.00.SAG.FK
	30 ●		XP-endo Rise Shaper	S1.XB0.00.0AV.FK	S1.XB0.00.0AS.FK	
31	15 O	0.04	XP-endo Rise Glider	-	S1.7B0.00.0HP.FK	3 x (XP-RG + XP-RS) S1.XB0.00.SAJ.FK
	30 ●		XP-endo Rise Shaper	S1.XB0.00.0AW.FK	S1.XB0.00.0AT.FK	



#### R-Motion®



NITI STERILE R

Length (mm)	ISO	Taper	Name	Blister 3 pcs
21	20 ●	0.04	R-Motion 20	S1.7B0.00.0HQ.FK
	25 ●	0.06	R-Motion 25	S1.7B0.00.0FR.FK
	30 ●	0.04	R-Motion 30	S1.7B0.00.0FI.FK
	40 ●		R-Motion 40	S1.7B0.00.0FJ.FK
	50 ●		R-Motion 50	S1.7B0.00.0FK.FK
25	20 ●	0.04	R-Motion 20	S1.7B0.00.0HR.FK
	25 ●	0.06	R-Motion 25	S1.7B0.00.0FS.FK
	30 ●	0.04	R-Motion 30	S1.7B0.00.0FL.FK
	40 ●		R-Motion 40	S1.7B0.00.0FM.FK
	50 ●		R-Motion 50	S1.7B0.00.0FN.FK
31	20 ●	0.04	R-Motion 20	S1.7B0.00.0HS.FK
	25 ●	0.06	R-Motion 25	S1.7B0.00.0FT.FK
	30 ●	0.04	R-Motion 30	S1.7B0.00.0FO.FK
	40 ●		R-Motion 40	S1.7B0.00.0FP.FK
	50 ●		R-Motion 50	S1.7B0.00.0FQ.FK



For specific R-Motion glide path files, see [page 9](#).



## RACE® EVO 4%



Length (mm)	ISO	Taper	Name	Blister 6 pcs	Sequence 3 pcs
21	15 ○	0.04	RE1	S1.7B0.00.0DG.FK	RE1 + RE2 + RE3 S1.7B0.00.SCV.FK
	25 ●		RE2	S1.7B0.00.0DI.FK	
	30 ●		RE3	S1.7B0.00.0DK.FK	
	40 ●		-	S1.7B0.00.0DM.FK	-
	50 ●		-	S1.7B0.00.0DO.FK	-
25	15 ○		RE1	S1.7B0.00.0FZ.FK	RE1 + RE2 + RE3 S1.7B0.00.SCW.FK
	25 ●		RE2	S1.7B0.00.0GB.FK	
	30 ●		RE3	S1.7B0.00.0GC.FK	
	40 ●		-	S1.7B0.00.0GE.FK	-
	50 ●		-	S1.7B0.00.0GF.FK	-
31	15 ○		RE1	S1.7B0.00.0DH.FK	RE1 + RE2 + RE3 S1.7B0.00.SCX.FK
	25 ●		RE2	S1.7B0.00.0DJ.FK	
	30 ●		RE3	S1.7B0.00.0DL.FK	
	40 ●		-	S1.7B0.00.0DN.FK	-
	50 ●		-	S1.7B0.00.0DP.FK	-



## RACE® EVO 6%



Length (mm)	ISO	Taper	Name	Blister 6 pcs	Sequence 3 pcs
21	15 ○	0.04	RE1	S1.7B0.00.0DG.FK	RE1 + RE2 + RE3 S1.7B0.00.SCY.FK
	25 ●		RE2	S1.7B0.00.0DI.FK	
	25 ●	0.06	RE3	S1.7B0.00.0DQ.FK	-
	30 ●		-	S1.7B0.00.0DV.FK	-
25	15 ○	0.04	RE1	S1.7B0.00.0FZ.FK	RE1 + RE2 + RE3 S1.7B0.00.SCZ.FK
	25 ●		RE2	S1.7B0.00.0GB.FK	
	25 ●	0.06	RE3	S1.7B0.00.0GH.FK	-
	30 ●		-	S1.7B0.00.0GI.FK	-
	35 ●		-	S1.7B0.00.0GJ.FK	-
31	15 ○	0.04	RE1	S1.7B0.00.0DH.FK	RE1 + RE2 + RE3 S1.7B0.00.SDA.FK
	25 ●		RE2	S1.7B0.00.0DJ.FK	
	25 ●	0.06	RE3	S1.7B0.00.0DU.FK	-
	30 ●		-	S1.7B0.00.0DW.FK	-
	35 ●		-	S1.7B0.00.0GU.FK	-




## RACE® EVO



Length (mm)	ISO	Taper	Blister 6 pcs	
21	30 ●	0.02	S1.7B0.00.0HB.FK	
	20 ●		S1.7B0.00.0HF.FK	
	25 ●		S1.7B0.00.0DI.FK	
	30 ●		S1.7B0.00.0DK.FK	
	35 ●		S1.7B0.00.0HH.FK	
	25	40 ●	0.04	S1.7B0.00.0DM.FK
		50 ●		S1.7B0.00.0DO.FK
		25 ●		S1.7B0.00.0DQ.FK
		30 ●		S1.7B0.00.0DV.FK
		35 ●		S1.7B0.00.0DX.FK
31		40 ●	0.06	S1.7B0.00.0HL.FK
		30 ●		S1.7B0.00.0FX.FK
		20 ●		S1.7B0.00.0GA.FK
		25 ●		S1.7B0.00.0GB.FK
		30 ●		S1.7B0.00.0GC.FK
	25	35 ●	0.04	S1.7B0.00.0GD.FK
		40 ●		S1.7B0.00.0GE.FK
		50 ●		S1.7B0.00.0GF.FK
		25 ●		S1.7B0.00.0GH.FK
		30 ●		S1.7B0.00.0GI.FK
31		35 ●	0.06	S1.7B0.00.0GJ.FK
		40 ●		S1.7B0.00.0GK.FK
		30 ●		S1.7B0.00.0HC.FK
		20 ●		S1.7B0.00.0HG.FK
		25 ●		S1.7B0.00.0DJ.FK
	25	30 ●	0.04	S1.7B0.00.0DL.FK
		35 ●		S1.7B0.00.0HI.FK
		40 ●		S1.7B0.00.0DN.FK
		50 ●		S1.7B0.00.0DP.FK
		25 ●		S1.7B0.00.0DU.FK
31		30 ●	0.06	S1.7B0.00.0DW.FK
		35 ●		S1.7B0.00.0GU.FK
		40 ●		S1.7B0.00.0HM.FK
		30 ●		S1.7B0.00.0DU.FK
		35 ●		S1.7B0.00.0DW.FK
	25	40 ●	0.06	S1.7B0.00.0GU.FK
		30 ●		S1.7B0.00.0DU.FK
		35 ●		S1.7B0.00.0DW.FK
		40 ●		S1.7B0.00.0HM.FK
		30 ●		S1.7B0.00.0DU.FK




For specific RACE EVO glide path files, see [page 9](#).



# B. Files

## 04

### Cleaning



Finisher



NI Ti STERILE R

Length (mm)	ISO	Taper	Blister 3 pcs
21	25	0.00	S1.XB0.00.0AX.FK
25			S1.XB0.00.0AY.FK
31			S1.XB0.00.0AZ.FK



Finisher R



NI Ti STERILE R

Length (mm)	ISO	Taper	Blister 3 pcs
21	30	0.00	S1.XB0.00.0BA.FK
25			S1.XB0.00.0BB.FK
31			S1.XB0.00.0BC.FK







# B. Files

# 05

## Retreatment



XP-endo® Rise  
Retreatment



NI Ti STERILE R

Length (mm)	ISO	Taper	Name	Sequence 3 pcs
21		0.10	DR1	DR1 + XP-RS + Finisher R S1.XB0.00.SAK.FK
		0.04	XP-endo Rise Shaper	
		0.00	Finisher R	
25	30	0.10	DR1	DR1 + XP-RS + Finisher R S1.XB0.00.SAH.FK
		0.04	XP-endo Rise Shaper	
		0.00	Finisher R	
31		0.10	DR1	DR1 + XP-RS + Finisher R S1.XB0.00.SAL.FK
		0.04	XP-endo Rise Shaper	
		0.00	Finisher R	



D-Race



NI Ti STERILE R

Length (mm)	ISO	Taper	Name	Blister 6 pcs	Sequences 4 pcs
15	30	0.10	DR1	S1.7B0.00.0FB.FK	2 × (DR1 + DR2) S1.7B0.00.SAI.FK
25	25	0.04	DR2	S1.7B0.00.0FA.FK	




# C. Hand Files



## 01

K-Files, H-Files, Reamers, F-Search, Barbed broaches




K-Files



Length (mm)	ISO	Taper	Blister 6 pcs	Assortment 6 pcs
21	6	0.02	S1.1B4.00.0AA.FK	ISO 6-25 S1.1B4.00.SAJ.FK  ISO 15-40 S1.1B4.00.SAA.FK  ISO 45-80 S1.1B4.00.SAD.FK
	8		S1.1B4.00.0AB.FK	
	10		S1.1B4.00.0AC.FK	
	15		S1.1B4.00.0AD.FK	
	20		S1.1B4.00.0AE.FK	
	25		S1.1B4.00.0AF.FK	
	30		S1.1B4.00.0AG.FK	
	35		S1.1B4.00.0AH.FK	
	40		S1.1B4.00.0AI.FK	
	45			
	50			
	55			
	60			
	70			
	80			
25	6	0.02	S1.1B4.00.0AJ.FK	ISO 6-25 S1.1B4.00.SAK.FK  ISO 08-15 S1.1B4.00.SAM.FK (2 x 3 pcs)  ISO 15-40 S1.1B4.00.SAB.FK  ISO 45-80 S1.1B4.00.SAF.FK
	8		S1.1B4.00.0AK.FK	
	10		S1.1B4.00.0AL.FK	
	15		S1.1B4.00.0AM.FK	
	20		S1.1B4.00.0AN.FK	
	25		S1.1B4.00.0AO.FK	
	30		S1.1B4.00.0AP.FK	
	35		S1.1B4.00.0AQ.FK	
	40		S1.1B4.00.0AR.FK	
	45			
	50			
	55			
	60			
	70			
	80			
31	6	0.02	S1.1B4.00.0AS.FK	ISO 6-25 S1.1B4.00.SAL.FK  ISO 15-40 S1.1B4.00.SAC.FK  ISO 45-80 S1.1B4.00.SAH.FK
	8		S1.1B4.00.0AT.FK	
	10		S1.1B4.00.0AU.FK	
	15		S1.1B4.00.0AV.FK	
	20		S1.1B4.00.0AW.FK	
	25		S1.1B4.00.0AX.FK	
	30		S1.1B4.00.0AY.FK	
	35		S1.1B4.00.0AZ.FK	
	40		S1.1B4.00.0BA.FK	
	45			
	50			
	55			
	60			
	70			
	80			




K-Files Flex



Length (mm)	ISO	Taper	Plastic Box 6 pcs SSt	Assortment 6 pcs SSt	Plastic Box 6 pcs NiTi	Assortment 6 pcs NiTi
21	8	0.02	01.554.21.008.CG	-	-	-
	10		01.554.21.010.CG	-	-	-
	15		01.554.21.015.CG	-	01.754.21.015.CG	-
	20		01.554.21.020.CG	-	01.754.21.020.CG	-
	25		01.554.21.025.CG	ISO 15-40 01.554.21.S15.CG	01.754.21.025.CG	ISO 15-40 01.754.21.S15.CG
	30		01.554.21.030.CG	-	01.754.21.030.CG	-
	35		01.554.21.035.CG	-	01.754.21.035.CG	-
	40		01.554.21.040.CG	-	01.754.21.040.CG	-
25	8	0.02	01.554.25.008.CG	-	-	-
	10		01.554.25.010.CG	-	-	-
	15		01.554.25.015.CG	-	01.754.25.015.CG	-
	20		01.554.25.020.CG	-	01.754.25.020.CG	-
	25		01.554.25.025.CG	ISO 15-40 01.554.25.S15.CG	01.754.25.025.CG	ISO 15-40 01.754.25.S15.CG
	30		01.554.25.030.CG	-	01.754.25.030.CG	-
	35		01.554.25.035.CG	-	01.754.25.035.CG	-
	40		01.554.25.040.CG	-	01.754.25.040.CG	-
31	8	0.02	01.554.31.008.CG	-	-	-
	10		01.554.31.010.CG	-	-	-
	15		01.554.31.015.CG	-	-	-
	20		01.554.31.020.CG	-	-	-
	25		01.554.31.025.CG	ISO 15-40 01.554.31.S15.CG	-	-
	30		01.554.31.030.CG	-	-	-
	35		01.554.31.035.CG	-	-	-
	40		01.554.31.040.CG	-	-	-



H-Files



Length (mm)	ISO	Taper	Plastic Box 6 pcs	Assortment 6 pcs		
21	8	0.02	01.153.21.008.CG	-		
	10		01.153.21.010.CG	-		
	15		01.153.21.015.CG	-		
	20		01.153.21.020.CG	-		
	25		01.153.21.025.CG	-		
	30		01.153.21.030.CG	-		
	35		01.153.21.035.CG	ISO 15-40 01.153.21.S15.CG		
	40		01.153.21.040.CG	-		
	45		01.153.21.045.CG	ISO 45-80 01.153.21.S45.CG		
	50		01.153.21.050.CG	-		
	55		01.153.21.055.CG	-		
	60		01.153.21.060.CG	-		
	70		01.153.21.070.CG	-		
	80		01.153.21.080.CG	-		
	25		8	0.02	01.153.25.008.CG	-
			10		01.153.25.010.CG	-
15		01.153.25.015.CG	-			
20		01.153.25.020.CG	-			
25		01.153.25.025.CG	-			
30		01.153.25.030.CG	-			
35		01.153.25.035.CG	ISO 15-40 01.153.25.S15.CG			
40		01.153.25.040.CG	-			
45		01.153.25.045.CG	ISO 45-80 01.153.25.S45.CG			
50		01.153.25.050.CG	-			
55		01.153.25.055.CG	-			
60		01.153.25.060.CG	-			
70		01.153.25.070.CG	-			
80		01.153.25.080.CG	-			
31		8	0.02		01.153.31.008.CG	-
		10			01.153.31.010.CG	-
	15	01.153.31.015.CG		-		
	20	01.153.31.020.CG		-		
	25	01.153.31.025.CG		-		
	30	01.153.31.030.CG		-		
	35	01.153.31.035.CG		ISO 15-40 01.153.31.S15.CG		
	40	01.153.31.040.CG		-		
	45	01.153.31.045.CG		ISO 45-80 01.153.31.S45.CG		
	50	01.153.31.050.CG		-		
	55	01.153.31.055.CG		-		
	60	01.153.31.060.CG		-		
	70	01.153.31.070.CG		-		
	80	01.153.31.080.CG		-		



Reamers



Length (mm)	ISO	Taper	Plastic Box 6 pcs	Assortment 6 pcs
21	6	0.02	01.160.21.006.CG	ISO 6-25 01.160.21.S06.CG  ISO 15-40 01.160.21.S15.CG  ISO 45-80 01.160.21.S45.CG
	8		01.160.21.008.CG	
	10		01.160.21.010.CG	
	15		01.160.21.015.CG	
	20		01.160.21.020.CG	
	25		01.160.21.025.CG	
	30		01.160.21.030.CG	
	35		01.160.21.035.CG	
	40		01.160.21.040.CG	
	45		01.160.21.045.CG	
	50		01.160.21.050.CG	
	55		01.160.21.055.CG	
	60		01.160.21.060.CG	
70	01.160.21.070.CG			
80	01.160.21.080.CG			
25	6	0.02	01.160.25.006.CG	ISO 6-25 01.160.25.S06.CG  ISO 15-40 01.160.25.S15.CG  ISO 45-80 01.160.25.S45.CG
	8		01.160.25.008.CG	
	10		01.160.25.010.CG	
	15		01.160.25.015.CG	
	20		01.160.25.020.CG	
	25		01.160.25.025.CG	
	30		01.160.25.030.CG	
	35		01.160.25.035.CG	
	40		01.160.25.040.CG	
	45		01.160.25.045.CG	
	50		01.160.25.050.CG	
	55		01.160.25.055.CG	
	60		01.160.25.060.CG	
70	01.160.25.070.CG			
80	01.160.25.080.CG			
31	6	0.02	01.160.31.006.CG	ISO 15-40 01.160.31.S15.CG  ISO 45-80 01.160.31.S45.CG
	8		01.160.31.008.CG	
	10		01.160.31.010.CG	
	15		01.160.31.015.CG	
	20		01.160.31.020.CG	
	25		01.160.31.025.CG	
	30		01.160.31.030.CG	
	35		01.160.31.035.CG	
	40		01.160.31.040.CG	
	45		01.160.31.045.CG	
	50		01.160.31.050.CG	
	55		01.160.31.055.CG	
	60		01.160.31.060.CG	
70	01.160.31.070.CG			
80	01.160.31.080.CG			



F-Search



Length (mm)	ISO	Taper	Plastic Box 6 pcs
21	20	0.05	03.176.21.020.CG



Barbed Broaches



Length (mm)	ISO	Taper	Plastic Box 144 pcs*	Assortment 144 pcs**
20	25	-	03.050.31.025.CT	ISO 25-60 03.050.31.S25.CT
	30		03.050.31.030.CT	
	35		03.050.31.035.CT	
	40		03.050.31.040.CT	
	50		03.050.31.050.CT	
	60		03.050.31.060.CT	



\* 12 boxes of 12 instruments  
\*\* 2 boxes of 12 instruments of each ISO reference




# D. Obturation



# 01

## Bioceramic obturation



TotalFill® Range /  
Gamme TotalFill®



Name	Content	Ref.
TotalFill® BC Sealer™	Syringe 2g + 15 TotalFill Tips	08.650.00.007.XX
TotalFill® BC Sealer™ HiFlow	Syringe 1.5g + 15 TotalFill Tips	08.650.00.008.XX
TotalFill® BC RRM™ Paste	Syringe 1g + 7 TotalFill Tips	08.650.00.003.XX
TotalFill® BC RRM™ Putty	Jar 2.5g	08.650.00.004.XX
TotalFill® BC RRM™ Fast Set Putty	Syringe 0.3g	08.650.00.006.XX
TotalFill® BC Tips™	15 Tips	08.650.00.002.XX
TotalFill® BC Pellets™	100 Pellets	08.650.00.005.XX




TotalFill® BC Points™ HiFlow



ISO	Taper	Ref.	Quantity	Assortment
15 ○	0.04 ■	08.650.04.015.HF	x 60 pcs	-
20 ●		08.650.04.020.HF		
25 ●		08.650.04.025.HF		
30 ●		08.650.04.030.HF		
35 ●		08.650.04.035.HF		
40 ●		08.650.04.040.HF		
45 ○		08.650.04.045.HF		
55 ●		08.650.04.055.HF		



TotalFill® BC Points™



ISO	Taper	Ref.	Quantity	Assortment
25 ●	0.02 ■	08.650.02.025.RE	x 100 pcs	ISO 15-60 08.650.02.S60.RB (x 300 pcs)
30 ●		08.650.02.030.RE		
35 ●		08.650.02.035.RE		
40 ●		08.650.02.040.RE		
45 ○		08.650.02.045.RE		
50 ●		08.650.02.050.RE		
55 ●		08.650.02.055.RE		
60 ●		08.650.02.060.RE		
15 ○	0.04 ■	08.650.04.015.RE	x 60 pcs	ISO 15-60 08.650.04.S60.RB (x 300 pcs)
20 ●		08.650.04.020.RE		
25 ●		08.650.04.025.RE		
30 ●		08.650.04.030.RE		
35 ●		08.650.04.035.RE		
40 ●		08.650.04.040.RE		
45 ○		08.650.04.045.RE		
50 ●		08.650.04.050.RE		
55 ●	08.650.04.055.RE			
60 ●	08.650.04.060.RE			
70 ●	08.650.04.070.RE			
80 ●	08.650.04.080.RE			
15 ○	0.06 ■	08.650.06.015.RE	x 60 pcs	ISO 15-50 08.650.06.S60.RB (x 240 pcs)
20 ●		08.650.06.020.RE		
25 ●		08.650.06.025.RE		
30 ●		08.650.06.030.RE		
35 ●		08.650.06.035.RE		
40 ●		08.650.06.040.RE		
45 ○	08.650.06.045.RE			
50 ●	08.650.06.050.RE			





BC Universal



Name	Content	Ref.
BC Universal Sealer	Syringe 2g + 20 Disposable tips + Syringe Holder	08.650.00.013.XX
BC Universal RRM	10 Capsules - 0.25g each	08.650.00.014.XX



02

## Paper points &amp; Gutta percha

Gutta Percha &amp; Paper Points




ISO	Taper	Ref. Gutta Percha	Ref. Paper Points <i>Sterile</i>	Quantity	Assortment Gutta Percha	Assortment Paper Point			
15 ○	0.02	08.600.02.015.RE	08.610.02.015.BC	x 120 pcs (Gutta Percha)  x 200 pcs (Paper Points)	ISO 15-40 08.600.02.S15.RE	ISO 15-40 08.610.02.S15.BC			
20 ●		08.600.02.020.RE	08.610.02.020.BC						
25 ●		08.600.02.025.RE	08.610.02.025.BC						
30 ●		08.600.02.030.RE	08.610.02.030.BC						
35 ●		08.600.02.035.RE	08.610.02.035.BC						
40 ●		08.600.02.040.RE	08.610.02.040.BC						
45 ○		08.600.02.045.RE	08.610.02.045.BC						
50 ●		08.600.02.050.RE	08.610.02.050.BC						
55 ●		08.600.02.055.RE	08.610.02.055.BC						
60 ●		08.600.02.060.RE	08.610.02.060.BC						
70 ●	08.600.02.070.RE	08.610.02.070.BC	x 60 pcs (Gutta Percha)  x 200 pcs (Paper Points)	ISO 15-40 08.600.04.S15.RE	ISO 15-40 08.610.04.S15.BC				
15 ○	08.600.04.015.RE	08.610.04.015.BC							
20 ●	08.600.04.020.RE	08.610.04.020.BC							
25 ●	08.600.04.025.RE	08.610.04.025.BC							
30 ●	08.600.04.030.RE	08.610.04.030.BC							
35 ●	08.600.04.035.RE	08.610.04.035.BC							
40 ●	08.600.04.040.RE	08.610.04.040.BC							
50 ●	08.600.04.050.RE	08.610.04.050.BC							
15 ○	0.06	08.600.06.015.RE				08.610.06.015.BC	x 60 pcs (Gutta Percha)  x 200 pcs (Paper Points)	ISO 15-40 08.600.06.S15.RE	ISO 15-40 08.610.06.S15.BC
20 ●		08.600.06.020.RE				08.610.06.020.BC			
25 ●		08.600.06.025.RE	08.610.06.025.BC						
30 ●		08.600.06.030.RE	08.610.06.030.BC						
35 ●		08.600.06.035.RE	08.610.06.035.BC						
40 ●		08.600.06.040.RE	08.610.06.040.BC						



# 03

## Paste fillers

Lentulo & Spring Paste Fillers



Length (mm)	ISO	Plastic Box Lentulo 4 pcs	Plastic Box Spring 4 pcs	Assortment Lentulo 4 pcs	Assortment Spring 4 pcs
21	25 ●	02.171.21.025.GM	02.190.21.025.GM	ISO 25-40 02.171.21.S25.GM	ISO 25-40 02.190.21.S25.GM
	30 ●	02.171.21.030.GM	02.190.21.030.GM		
	35 ●	02.171.21.035.GM	02.190.21.035.GM		
	40 ●	02.171.21.040.GM	02.190.21.040.GM		
25	25 ●	02.171.25.025.GM	02.190.25.025.GM	ISO 25-40 02.171.25.S25.GM	ISO 25-40 02.190.25.S25.GM
	30 ●	02.171.25.030.GM	02.190.25.030.GM		
	35 ●	02.171.25.035.GM	02.190.25.035.GM		
	40 ●	02.171.25.040.GM	02.190.25.040.GM		
29	25 ●	02.171.29.025.GM	02.190.29.025.GM	ISO 25-40 02.171.29.S25.GM	ISO 25-40 02.190.29.S25.GM
	30 ●	02.171.29.030.GM	02.190.29.030.GM		
	35 ●	02.171.29.035.GM	02.190.29.035.GM		
	40 ●	02.171.29.040.GM	02.190.29.040.GM		



# 04

## Finger Spreaders & Pluggers



Spreaders



Length (mm)	ISO	Taper	Plastic Box 6 pcs SSt	Assortment 6 pcs SSt	Plastic Box 6 pcs NiTi	Assortment 6 pcs NiTi
21	15 ○	0.02	02.111.21.015.CG	ISO 15-40 02.111.21.S15.CG	02.311.21.015.CG	ISO 15-40 02.311.21.S15.CG
	20 ●		02.111.21.020.CG		02.311.21.020.CG	
	25 ●		02.111.21.025.CG		02.311.21.025.CG	
	30 ●		02.111.21.030.CG		02.311.21.030.CG	
	35 ●		02.111.21.035.CG		02.311.21.035.CG	
	40 ●		02.111.21.040.CG		02.311.21.040.CG	
25	20 ●	0.04	-	-	02.312.21.020.CG	-
	15 ○	0.02	02.111.25.015.CG	ISO 15-40 02.111.25.S15.CG	02.311.25.015.CG	ISO 15-40 02.311.25.S15.CG
	20 ●		02.111.25.020.CG		02.311.25.020.CG	
	25 ●		02.111.25.025.CG		02.311.25.025.CG	
	30 ●		02.111.25.030.CG		02.311.25.030.CG	
	35 ●		02.111.25.035.CG		02.311.25.035.CG	
40 ●	02.111.25.040.CG		02.311.25.040.CG			



Sensipast



Length (mm)	ISO	Plastic Box 4 pcs	Assortment 4 pcs
21	25 ●	02.180.21.025.GS	ISO 25-40 02.180.21.S25.GS
	30 ●	02.180.21.030.GS	
	35 ●	02.180.21.035.GS	
	40 ●	02.180.21.040.GS	
25	25 ●	02.180.25.025.GS	ISO 25-40 02.180.25.S25.GS
	30 ●	02.180.25.030.GS	
	35 ●	02.180.25.035.GS	
	40 ●	02.180.25.040.GS	
29	25 ●	02.180.29.025.GS	ISO 25-40 02.180.29.S25.GS
	30 ●	02.180.29.030.GS	
	35 ●	02.180.29.035.GS	
	40 ●	02.180.29.040.GS	




Pluggers



Length (mm)	ISO	Taper	Plastic Box 6 pcs SSt	Assortment 6 pcs SSt	Plastic Box 6 pcs NiTi	Assortment 6 pcs NiTi
21	15 ○	0.02	02.110.21.015.CG	ISO 15-40 02.110.21.S15.CG	02.310.21.015.CG	ISO 15-40 02.310.21.S15.CG
	20 ●		02.110.21.020.CG		02.310.21.020.CG	
	25 ●		02.110.21.025.CG		02.310.21.025.CG	
	30 ●		02.110.21.030.CG		02.310.21.030.CG	
	35 ●		02.110.21.035.CG		02.310.21.035.CG	
	40 ●		02.110.21.040.CG		02.310.21.040.CG	
25	15 ○	0.02	02.110.25.015.CG	ISO 15-40 02.110.25.S15.CG	02.310.25.015.CG	ISO 15-40 02.310.25.S15.CG
	20 ●		02.110.25.020.CG		02.310.25.020.CG	
	25 ●		02.110.25.025.CG		02.310.25.025.CG	
	30 ●		02.110.25.030.CG		02.310.25.030.CG	
	35 ●		02.110.25.035.CG		02.310.25.035.CG	
	40 ●		02.110.25.040.CG		02.310.25.040.CG	



# E. Accessories



## SMD



Color	Quantity	Perforation (ø, mm)	Ref.
Transparent	12	No	06.030.03.010.TR
		0.5	06.030.05.010.TR
		0.9	06.030.09.010.TR
	100	No	06.030.03.100.TR
		0.5	06.030.05.100.TR
		0.9	06.030.09.100.TR



## SMD dispensers



Content	Color	Perforation (ø, mm)	Ref.
1 Dispenser + 100 SMD	Transparent	0.5	06.031.05.000.TR



## Endo stops



Color	Quantity	Perforation	Ref.
Yellow	12	No	06.010.00.010.JA
Red			06.010.00.010.RO
Green			06.010.00.010.VE
Black			06.010.00.010.NO
Yellow	100	No	06.010.00.100.JA
Red			06.010.00.100.RO
Green			06.010.00.100.VE
Black			06.010.00.100.NO
Red	12	Yes	06.010.01.010.RO
Blue			06.010.01.010.BL
Green			06.010.01.010.VE
Black			06.010.01.010.NO
Red	100	Yes	06.010.01.100.RO
Blue			06.010.01.100.BL
Green			06.010.01.100.VE
Black			06.010.01.100.NO




## Endo stands



Name	Color	Arch (mm)	Ref.
Endo stand Freestyle	Orange	12	08.096.01.000.OR



## Pins



Color	Quantity	Dimension (mm)	Ref.
Transparent	100	1	05.300.00.101.TR
	50	2	05.300.00.052.TR



# F. Restoration




## Matrixes



Thickness (mm)	Quantity	N°	Ref. Premolar	N°	Ref. Molar
0.03	12	208	04.208.00.013.AC	206	04.206.00.013.AC
		209	04.209.00.013.AC	207	04.207.00.013.AC
		210	04.210.00.013.AC	309	04.309.00.013.AC
0.05		208	04.208.00.015.AC	206	04.206.00.015.AC
		209	04.209.00.015.AC	207	04.207.00.015.AC
		210	04.210.00.015.AC	309	04.309.00.015.AC
0.03	100	208	04.208.00.103.AC	206	04.206.00.103.AC
		209	04.209.00.103.AC	207	04.207.00.103.AC
		210	04.210.00.103.AC	309	04.309.00.103.AC
0.05		208	04.208.00.105.AC	206	04.206.00.105.AC
		209	04.209.00.105.AC	207	04.207.00.105.AC
		210	04.210.00.105.AC	309	04.309.00.105.AC



## Matrixes



Thickness (mm)	Quantity	N°	Ref. Premolar	N°	Ref. Molar
0.05	12	203	04.203.00.015.AC	200	04.200.00.015.AC
		204	04.204.00.015.AC	201	04.201.00.015.AC
		205	04.205.00.015.AC	202	04.202.00.015.AC
	100	203	04.203.00.105.AC	-	-
		204	04.204.00.105.AC	201	04.201.00.105.AC
		205	04.205.00.105.AC	202	04.202.00.105.AC










## Matrix Bands



Length (m)	Thickness (mm)	Width (mm)	Ref. 1 pcs
1	0.03	5	04.059.01.503.AC
		6	04.059.01.603.AC
		7	04.059.01.703.AC
	0.05	5	04.059.01.505.AC
		6	04.059.01.605.AC
		7	04.059.01.705.AC
5	0.03	5	04.059.05.503.AC
		6	04.059.05.603.AC
		7	04.059.05.703.AC
	0.05	5	04.059.05.505.AC
		6	04.059.05.605.AC
		7	04.059.05.705.AC
10	0.03	5	04.059.10.503.AC
		6	04.059.10.603.AC
		7	04.059.10.703.AC
	0.05	5	04.059.10.505.AC
		6	04.059.10.605.AC
		7	04.059.10.705.AC



## G. Legend

Symbol	Description
	Continuous rotation
	Reciprocating motion (counterclockwise)
	Hand held instrument
	Sterilized by irradiation
	Non Sterile
	Nickel titanium
	Stainless steel
	Silicone
	Aluminium
	Plastic
	More information (FKG Website)
	Protocol card
	Tutorial video (YouTube)

Device name	CE marking	Manufacturer
Barbed Broaches D-Race Finisher Finisher R Gates Lentulo Peeso PreRace Race	RACE® EVO R-Motion® Sensipast Spring Paste Fillers XP-endo® Rise Glider XP-endo® Rise Sequence XP-endo® Rise Shaper XP-endo® Rise Retreatment XS Gates	CE 1639
H-Files K-Files	K-Files Flex Reamers	CE/CE 1639
Calcinable pins Endo stands Endo stops F-Search Matrix bands	Matrixes Pluggers SMD SMD dispensers Spreaders	CE
Rooter® X3000		CE 0197
Rooter® Universal		CE 0197
TotalFill® BC Sealer™ TotalFill® BC Sealer™ HiFlow TotalFill® BC RRM™ Paste TotalFill® BC RRM™ Putty TotalFill® BC RRM™ Fast Set Putty TotalFill® BC Tips™		CE 0413
TotalFill® BC Pellets™ TotalFill® BC Points™ HiFlow TotalFill® BC Points™		CE 0197
BC Universal Sealer BC Universal RRM		CE 1639
Gutta Percha Paper Points		CE 1639

**FKG Dentaire Sàrl**  
Le Crêt-du-Loche 4  
2322 Le Crêt-du-Loche  
Switzerland

**Guilin Woodpecker Medical Instrument Co., Ltd.**  
Information Industrial Park,  
Guilin National High-Tech Zone,  
Guilin, 541004 Guangxi, P.R.  
China


**Foshan COXO Medical Instrument Co., Ltd.**  
No. 17, Guangming Ave.,  
New Light Source Industrial  
Base, Nanhai National High-  
tech Zone, Foshan 528226,  
Guangdong P.R. China

**Innovative BioCeramix Inc.**  
101-8218 North Fraser Way,  
Burnaby, British Columbia V3N  
0E9, Canada

**Sure Dent Corporation**  
# 809, 52, Sagimakgol-ro,  
Jungwon-gu, Seongnam-si,  
Gyeonggi-do, 13210, Republic  
of Korea

**Vericom Co., Ltd.**  
48, Toegyegongdan 1-gil,  
Chuncheon-Si, Gangwon-Do,  
24427, Republic of Korea

**Meta Biomed Co., Ltd.**  
270, Osongsaengmyeong 1-ro,  
Osong-eup, Heungdeok-gu,  
Cheongju-si, Chungcheongbuk-  
do, Republic of Korea





**FKG Dentaire Sàrl**  
Le Crêt-du-Loche 4  
2322 Le Crêt-du-Loche  
Switzerland  
T +41 32 924 22 44

[info@fkg.ch](mailto:info@fkg.ch)  
[www.fkg.ch](http://www.fkg.ch)



FKG Dentaire Sàrl  
Le Crêt-du-Loche 4  
2322 Le Crêt-du-Loche  
Switzerland

FKG REF. 99.000.00.07A.EN\_REV0\_2024-09-09 - © FKG Dentaire Sàrl - Pictures of our products are not contractual.